

Product datasheet for RC227802L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

alpha smooth muscle Actin (ACTA2) (NM 001141945) Human Tagged ORF Clone Lentiviral **Particle**

Product data:

Product Type: Lentiviral Particles

Product Name: alpha smooth muscle Actin (ACTA2) (NM_001141945) Human Tagged ORF Clone Lentiviral

Particle

Symbol: alpha smooth muscle Actin

Synonyms: ACTSA

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

mGFP Tag:

ACCN: NM 001141945

ORF Size: 1131 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(RC227802).

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional

amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA.

Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence

verification at a reduced cost. Please contact our customer care team at

custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 001141945.1





alpha smooth muscle Actin (ACTA2) (NM_001141945) Human Tagged ORF Clone Lentiviral Particle – RC227802L4V

RefSeq Size:1798 bpRefSeq ORF:1134 bp

Locus ID: 59

 UniProt ID:
 P62736

 Cytogenetics:
 10q23.31

Protein Pathways: Vascular smooth muscle contraction


MW: 42 kDa

Gene Summary: This gene encodes one of six different actin proteins. Actins are highly conserved proteins

that are involved in cell motility, structure, integrity, and intercellular signaling. The encoded protein is a smooth muscle actin that is involved in vascular contractility and blood pressure homeostasis. Mutations in this gene cause a variety of vascular diseases, such as thoracic aortic disease, coronary artery disease, stroke, and Moyamoya disease, as well as

multisystemic smooth muscle dysfunction syndrome. [provided by RefSeq, Sep 2017]

Product images:



[RC227802L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC227802L4V particle to overexpress human ACTA2-mGFP fusion protein.