

Product datasheet for **RC227808A1V**

Human Transcription Factor SP9 (SP9) (NM_001145250) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Transcription Factor SP9 (SP9) (NM_001145250) AAV Particle
Tag:	Myc-DDK
Symbol:	Transcription Factor SP9
Synonyms:	ZNF990
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

ORF Nucleotide
Sequence:

>RC227808 representing NM_001145250
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC227808 representing NM_001145250
Red=Cloning site Green=Tags(s)

MATSI LGEEPRFGTTPLAMLAATCNKIGNTSPLTTLPESSAFAGGFGHPWKRSSSSCNLGSSLSGFAVAT
GGRGSGGLAGGSGAANSFAFLASTSPTSSAFSSDYGLFNSAAAAAAGVSPQEAGGQSAFISKVHT
AADGLYPRVGMHPYESWYKSGFHSTLAAGEVTNGAASSWWDVHSSPGSWLEVNQNPAGGLQSSLHSGAPQ
ASLHSQLGTYNPDFSSLTHSAFSSSTGLGSSAAAAASHLLSTSQHLLAQDGFKPVLPSYSDSSAAVAAAAAS
AMISGAAAAAGGSSSARSARYSGRATCDCPNCQEAERLGPAGASLRKGLHSHCHIPGCGKYVGKTSHLK
AHLRWHTGERPFCVNLWFGCKRFRTRSDCLQRHLRTHTEKRFACPVCNKRFMRSDHLKSHKITHNGGGGG
KKGSSDSDTASNLETPRESPLDILHDSGVSAARAAAAAAGGKEASGPNDS

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Species:

Human

Serotype:

AAV-2

ACCN:

NM 001145250

ORF Size:

1452 bp

Buffer:

PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001145250.1</u>
RefSeq ORF:	1455 bp
Locus ID:	100131390
UniProt ID:	<u>P0CG40</u>
Cytogenetics:	2q31.1
MW:	48.7 kDa