

Product datasheet for **RC227825A1V**

Human ALS2CR1 (NIF3L1) (NM_001142355) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ALS2CR1 (NIF3L1) (NM_001142355) AAV Particle
Tag:	Myc-DDK
Symbol:	ALS2CR1
Synonyms:	ALS2CR1; CALS-7; MDS015
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC227825 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTGTCATCTTGCGTACGCCAGTCCCCACGACAGTCCGGTTGTAGATTCCCTGATCTGCAATTCTT
CCCGTTCTTCATGGATTTGAAGGCTCTCCTTTCTTCTTGAATGACTTTGCATCCCTCTCGTTTGCTGA
GAGTTGGGACAATGTTGGATTACTGGTGAACCAAGCCACCACATACTGTAAATACACTCTTCTGACC
AATGACCTGACTGAGGAAGTGATGGAGGAGGTGCTGCAAAGAAGGCAGACCTCATTCTCTCCTACCATC
CGCCTATCTTCCGACCCATGAAGCGCATAACCTGGAACACATGGAAGGAGCGCCTGGTGATCCGGGCTCT
GGAGAACAGAGTCGGTATCTACTCTCCTCATACAGCCTATGATGCTGCGCCCCAGGGCGTCAACAACTGG
TTGGCTAAAGGGCTTGAGCTTTGACCTCCAGGCCATACATCCTTCCAAAGCTCCCAACTACCTTACAG
AGGGAAACCACCGAGTAGAATTCAACGTAACTACACCAAGACCTGGACAAAGTCATGTCTGCAGTGAA
AGGAATTGACGGTGTTTCTGTCACTTCTTTTCTGCTAGGACTGGTAATGAGGAACAAACACGGATTAAT
CTGAATTGTAAGGCTTTGATGCAGGTGGTAGATTTTCTTCCCGGAACAAACAACCTTTATCAGA
AGACGGAAATTCTGTCACTGGAGAAGCCTTTGCTTCTACATACTGGAATGGGACGGTTATGCACACTGGA
TGAATCTGTCTCCCTGGCAACCATGATTGATCGAATAAAAAGACACCTAAACTATCTCATATTCGCTTA
GCCCTTGGGGTGGGAGAACCTTAGAGTCTCAAGTCAAAGTCGTGGCCCTGTGTGCTGGTTCTGGGAGCA
GCGTTCTGCAGGTGTTGAGGCTGACCTTTACCTCACAGGTGAGATGTCCCATCATGATACTTTGGATGC
TGCTTCCCAAGGAATAAATGTCATCCTCTGTGAACACAGCAACACTGAACGAGGCTTTCTTTCTGACCTT
CGAGATATGCTGGATTCTCACTTGGAGAATAAGATAAATTATCCTATCAGAGACTGACAGGGACCCCTC
TTCAGGTGGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC227825 protein sequence Red=Cloning site Green=Tags(s) MLSSCVRPVPTTVRFVDSLICNSSRSFMDLKALLSSLNDFASLSFAESWDNVGLLVEPSPHTVNTLFLT NDL TEEVMEEVLQKKADLILSYHPIIFRPMKRITWNTWKERLVIRALENRVGIYSPHTAYDAAPQGVNNW LAKGLGACTSRPIHPSKAPNYPTTEGHRVEFNVNYTQDLQVMSAVKGIDGVSVTSFSARTGNEEQTRIN LNCTQKALMQVVDLFLSRNKQLYQKTEILSLEKPLLLHTGMGRCLTLDSEVSLATMIDRIKRHLKLSHRL ALGVGRTLESQVKVVALCAGSGSSVLQGVADLYLTGEMSHHDTLDAASQGINVILCEHSNTERGFLSDL RDMLDSHLENKINIILSETDRDPLQVV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001142355
ORF Size:	1131 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001142355.1
RefSeq Size:	1446 bp
RefSeq ORF:	1053 bp
Locus ID:	60491
UniProt ID:	Q9GZT8
Cytogenetics:	2q33.1
MW:	42 kDa