

Product datasheet for **RC227829A1V**

Human P311 (NREP) (NM_001142476) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human P311 (NREP) (NM_001142476) AAV Particle
Tag: Myc-DDK
Symbol: P311
Synonyms: C5orf13; D4S114; P311; PRO1873; PTZ17; SEZ17
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC227829 representing NM_001142476
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTTATTACCCAGAACTCTTTGTCTGGGTCAGTCAAGAACCATTTCCAAACAAGGACATGGAGGGAA
GGCTTCCTAAGGGAAGACTTCCTGTCCCAAAGGAAGTGAACCGCAAGAAGACGATGAGACAAACGCTGC
CTCCCTGACTCCACTGGCAGCAGTGAAGTCCGCTCCCAAGAATCAGTTACCTCCACTTTTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC227829 representing NM_001142476
Red=Cloning site **Green**=Tags(s)

MVYYPELFWVVSQEPFPNKDMEGRLPKGRLPVPKEVNRKKNDETNAASLTPLGSELRSRPRISYLHFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001142476
ORF Size: 204 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001142476.1](#), [NP_001135948.1](#)

RefSeq Size: 2201 bp

RefSeq ORF: 207 bp

Locus ID: 9315

UniProt ID: [Q16612](#)

Cytogenetics: 5q22.1

MW: 7.9 kDa