

## Product datasheet for **RC227829L3V**

### **P311 (NREP) (NM\_001142476) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	P311 (NREP) (NM_001142476) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NREP
Synonyms:	C5orf13; D4S114; P311; PRO1873; PTZ17; SEZ17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001142476
ORF Size:	204 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227829).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001142476.1</a> , <a href="#">NP_001135948.1</a>
RefSeq Size:	2201 bp
RefSeq ORF:	207 bp
Locus ID:	9315
UniProt ID:	<a href="#">Q16612</a>
Cytogenetics:	5q22.1
MW:	7.9 kDa



[View online »](#)

**Gene Summary:**

May have roles in neural function. Ectopic expression augments motility of gliomas. Promotes also axonal regeneration (By similarity). May also have functions in cellular differentiation (By similarity). Induces differentiation of fibroblast into myofibroblast and myofibroblast ameboid migration. Increases retinoic-acid regulation of lipid-droplet biogenesis (By similarity). Down-regulates the expression of TGFB1 and TGFB2 but not of TGFB3 (By similarity). May play a role in the regulation of alveolar generation.[UniProtKB/Swiss-Prot Function]