

## Product datasheet for **RC227839L4V**

### Apoptosis inhibitor 5 (API5) (NM\_001142930) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoptosis inhibitor 5 (API5) (NM_001142930) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Apoptosis inhibitor 5
Synonyms:	AAC-11; AAC11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142930
ORF Size:	1572 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227839).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001142930.1</a>
RefSeq Size:	3782 bp
RefSeq ORF:	1575 bp
Locus ID:	8539
UniProt ID:	<a href="#">Q9BZZ5</a>
Cytogenetics:	11p12
Protein Families:	Druggable Genome
MW:	59 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an apoptosis inhibitory protein whose expression prevents apoptosis after growth factor deprivation. This protein suppresses the transcription factor E2F1-induced apoptosis and also interacts with, and negatively regulates Acinus, a nuclear factor involved in apoptotic DNA fragmentation. Its depletion enhances the cytotoxic action of the chemotherapeutic drugs. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]