

Product datasheet for **RC227850L4V**

ALS2CR1 (NIF3L1) (NM_001142356) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ALS2CR1
Synonyms:	ALS2CR1; CALS-7; MDS015
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142356
ORF Size:	855 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC227850).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001142356.1](#), [NP_001135828.1](#)

RefSeq ORF: 858 bp

Locus ID: 60491

UniProt ID: [Q9GZT8](#)

Cytogenetics: 2q33.1

MW: 32 kDa



Gene Summary:

May function as a transcriptional corepressor through its interaction with COPS2, negatively regulating the expression of genes involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]