

## Product datasheet for **RC227852L3V**

### **SKIP (INPP5K) (NM\_001135642) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SKIP (INPP5K) (NM_001135642) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SKIP
Synonyms:	MDCC4ID; PPS; SKIP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135642
ORF Size:	1347 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227852).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001135642.1</a>
RefSeq Size:	3391 bp
RefSeq ORF:	1119 bp
Locus ID:	51763
UniProt ID:	<a href="#">Q9BT40</a>
Cytogenetics:	17p13.3
Protein Families:	Druggable Genome, Phosphatase


[View online »](#)

<b>Protein Pathways:</b>	Inositol phosphate metabolism, Insulin signaling pathway, Metabolic pathways, Phosphatidylinositol signaling system
<b>MW:</b>	51.1 kDa
<b>Gene Summary:</b>	This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]