

Product datasheet for **RC227888A1V**

Human SSPN (NM_001135823) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SSPN (NM_001135823) AAV Particle
Tag:	Myc-DDK
Symbol:	SSPN
Synonyms:	DAGA5; KRAG; NSPN; SPN1; SPN2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC227888 representing NM_001135823 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCAAGAACAAGCAGCCACGCGCCAGCAGAGGCAGGGGGGCCCGCGCCGCGGACGCCGCTGGGC
CCGACGACATGGAGCCGAAGAAGGGCACGGGGGCCCCAAGGAGTGCGGGGAGGAGAGCCCCGGACCTG
CTGCGGCTGCCGTTCCCGCTGCTGCTCGCCCTGCTGCAGCTGGCCCTGGGCATCGCCGTGACCGTGGTG
GGCTTCCTCATGGCGAGCATCAGCTCCTCCCTGCTAGTCAGGGACACTCCATTTGGGCTGGGATCATTG
TCTGCTTAGTGGCCTATCTTGGCTTGTATGCTTTGTGTCTCATATCAGGTTGACGAACGGACATGTAT
TCAATTTTCTATGAAACTGTTATACTTTCTGCTGAGTGCCCTGGGCCTGACGGTCTGTGTCTGGCCGTG
GCCTTTGCCGCCACCACTATTTCGAGCTCACACAGTTTACCTGTGAGACCACACTCGACTCTTGCCAGT
GCAAACTGCCCTCCTCGGAGCCGCTCAGCAGGACCTTTGTTTACCGGGATGTGACGGACTGTACCAGCGT
CACTGGCACTTTCAAAGTGTCTTACTCATCCAGATGATTCTTAATTTGGTCTGCGGCCTTGTGTGCTTG
TTGGCTGCTTTGTGATGTGAAACATAGGTACCAGGTCTTCTATGTGGGTGTCAGGATATGCTCCCTCA
CGGCTTCCGAAGGCCCCAGCAAAGATC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC227888 representing NM_001135823 Red=Cloning site Green=Tags(s) MGKNKQPRGQQRQGGPPAADAAGPDDMEPKGTGAPKECGEEEPRTCCGCRFPLLLALLQLALGIAVTVV GFLMASISSLLVRDTPFWAGIIVCLVAYLGLFMLCVSYQVDERTCIQFSMKLLYFLLSALGLTVCVLAV AFAAHHYSQLTQFTCETTLDSQCCKLPSEPLSRTFVYRDVTDCTSVTGTFKLFLLIQMILNLVCGLVCL LACFVMWKHRYQVFYVGVRICSLTASEGPQQKI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001135823
ORF Size:	729 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001135823.1</u>
RefSeq ORF:	423 bp
Locus ID:	8082
UniProt ID:	<u>Q14714</u>
Cytogenetics:	12p12.1
MW:	26.62 kDa