

Product datasheet for **RC227927L4V**

CD133 (PROM1) (NM_001145851) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CD133 (PROM1) (NM_001145851) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PROM1 |
| Synonyms: | AC133; CD133; CORD12; MCDR2; MSTP061; PROML1; RP41; STGD4 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001145851 |
| ORF Size: | 2499 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC227927). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001145851.1 |
| RefSeq ORF: | 2502 bp |
| Locus ID: | 8842 |
| UniProt ID: | O43490 |
| Cytogenetics: | 4p15.32 |
| Protein Families: | Druggable Genome, ES Cell Differentiation/IPS, Transmembrane |
| MW: | 93.5 kDa |



[View online »](#)

Gene Summary:

This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation. Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]