

Product datasheet for **RC227937L4V**

FAIM3 (FCMR) (NM_001142472) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAIM3 (FCMR) (NM_001142472) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAIM3
Synonyms:	Fas apoptotic inhibitory molecule 3; OTTHUMP00000034619; regulator of Fas-induced apoptosis; TOSO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142472
ORF Size:	1170 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227937).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001142472.1 , NP_001135944.1
RefSeq Size:	3067 bp
RefSeq ORF:	1172 bp
Locus ID:	9214
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome, Transmembrane
MW:	43.1 kDa



[View online »](#)

Gene Summary:

Fc receptors specifically bind to the Fc region of immunoglobulins (Igs) to mediate the unique functions of each Ig class. FAIM3 encodes an Fc receptor for IgM (see MIM 147020) (Kubagawa et al., 2009 [PubMed 19858324]; Shima et al., 2010 [PubMed 20042454]).[supplied by OMIM, Jul 2010]