

## Product datasheet for **RC227982**

### LSM5 (NM\_001139499) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** LSM5 (NM\_001139499) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** LSM5  
**Synonyms:** YER146W  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC227982 representing NM\_001139499  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

**ATGGCGGCTAACGCTACTACCAACCCGTCGAGCTGCTGCCGTTAGAGCTTGTGGACAAATGTATAGGAT  
CAAGAATTCACATCGTGATGAAGAGTGATAAGGAAATTGTTGGTACTCTCTAGGATTTGATGACTTTGT  
CAATATGGTACTGGAAGATGTCACTGAGTTTGAATCACACCAGAAGGAAGAAGGATTACTAAATTAGAT  
CAGATTTTGCTAAATGGAATAATAACAATGCTGGTTCTGGAGGAGAAGGACCTGAAGTG**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA**

**Protein Sequence:** >RC227982 representing NM\_001139499  
**Red**=Cloning site **Green**=Tags(s)

MAANATTNPSQLLPLELVKDCIGSRIHIVMKSKEIVGTLGFDFFVNMVLEDVTEFEITPEGRRITKLD  
QILLNNGNITMLVPGGEGPEV

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Restriction Sites:** Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001139499

ORF Size: 276 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001139499.1</a></u> , <u><a href="#">NP_001132971.1</a></u>
<b>RefSeq Size:</b>	2451 bp
<b>RefSeq ORF:</b>	189 bp
<b>Locus ID:</b>	23658
<b>UniProt ID:</b>	<u><a href="#">Q9Y4Y9</a></u>
<b>Cytogenetics:</b>	7p14.3
<b>Protein Families:</b>	Stem cell - Pluripotency
<b>Protein Pathways:</b>	RNA degradation, Spliceosome
<b>MW:</b>	9.9 kDa
<b>Gene Summary:</b>	Sm-like proteins were identified in a variety of organisms based on sequence homology with the Sm protein family (see SNRPD2; MIM 601061). Sm-like proteins contain the Sm sequence motif, which consists of 2 regions separated by a linker of variable length that folds as a loop. The Sm-like proteins are thought to form a stable heteromer present in tri-snRNP particles, which are important for pre-mRNA splicing.[supplied by OMIM, Apr 2004]