

## Product datasheet for **RC227995L4V**

### CCDC85C (NM\_001144995) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC85C (NM_001144995) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC85C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001144995
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227995).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001144995.1</a>
RefSeq ORF:	1260 bp
Locus ID:	317762
UniProt ID:	<a href="#">A6NKD9</a>
Cytogenetics:	14q32.2
MW:	45 kDa
Gene Summary:	May play an important role in cortical development, especially in the maintenance of radial glia.[UniProtKB/Swiss-Prot Function]


[View online »](#)