

Product datasheet for **RC228023L3V**

Small integral membrane protein 6 (SMIM6) (NM_001162997) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Small integral membrane protein 6 (SMIM6) (NM_001162997) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Small integral membrane protein 6
Synonyms:	C17orf110; ELN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001162997
ORF Size:	186 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228023).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162997.1
RefSeq ORF:	189 bp
Locus ID:	100130933
UniProt ID:	P0DI80
Cytogenetics:	17q25.1
MW:	6.8 kDa


[View online »](#)