

Product datasheet for **RC228030A1V**

Human ASXL1 (NM_001164603) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human ASXL1 (NM_001164603) AAV Particle

Tag: Myc-DDK

Symbol: ASXL1

Synonyms: BOPS; MDS

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC228030 representing NM_001164603
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGGACAAACAGAAGAAGAAGAAGGAGCGCACGTGGGCCGAGGCCGCGCCTGGTATTAGAAAAC
ACTCGGATGCTCCAATGACACCAAAACAGATTCTGCAGGTCATAGAGGCAGAAGGACTAAAGGAAATGAG
AAGTGGACTTCCCCTCTCGCATGCCTCAATGCTATGCTACATTCCAATTCAAGAGGAGGAGGGGTTG
TTTTATAAACTGCCTGGCCGAATCAGCCTTTTCACGCTCAAGGTG

ACGCGTACGCGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228030 representing NM_001164603
Red=Cloning site Green=Tags(s)

MKDKQKKKKERTWAEAARLVLENYSAPMTPKQILQVIEAEGLEKEMRSGTSPLACLNAMLHSNSRGGGL
FYKLPGRISLFTLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001164603

ORF Size: 255 bp

Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001164603.1</u> , <u>NP_001158075.1</u>
RefSeq ORF:	258 bp
Locus ID:	171023
Cytogenetics:	20q11.21
MW:	9.4 kDa