

Product datasheet for **RC228030L2V**

ASXL1 (NM_001164603) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ASXL1
Synonyms:	BOPS; MDS
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001164603
ORF Size:	255 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228030).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001164603.1 , NP_001158075.1
RefSeq ORF:	258 bp
Locus ID:	171023
Cytogenetics:	20q11.21
MW:	9.4 kDa



Gene Summary:

This gene is similar to the Drosophila additional sex combs gene, which encodes a chromatin-binding protein required for normal determination of segment identity in the developing embryo. The protein is a member of the Polycomb group of proteins, which are necessary for the maintenance of stable repression of homeotic and other loci. The protein is thought to disrupt chromatin in localized areas, enhancing transcription of certain genes while repressing the transcription of other genes. The protein encoded by this gene functions as a ligand-dependent co-activator for retinoic acid receptor in cooperation with nuclear receptor coactivator 1. Mutations in this gene are associated with myelodysplastic syndromes and chronic myelomonocytic leukemia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]