

## Product datasheet for **RC228038L3V**

### **BANF2 (NM\_001159495) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BANF2 (NM_001159495) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BANF2
Synonyms:	BAF-L; BAF2; BAFL; C20orf179
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001159495
ORF Size:	291 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228038).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001159495.1</a> , <a href="#">NP_001152967.1</a>
RefSeq ORF:	294 bp
Locus ID:	140836
UniProt ID:	<a href="#">Q9H503</a>
Cytogenetics:	20p12.1
MW:	10.9 kDa
Gene Summary:	May play a role in BANF1 regulation and influence tissue-specific roles of BANF1. [UniProtKB/Swiss-Prot Function]



[View online »](#)