

Product datasheet for **RC228039L4V**

BRI3 (NM_001159491) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	BRI3
Synonyms:	I3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159491
ORF Size:	294 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228039).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001159491.1
RefSeq Size:	699 bp
RefSeq ORF:	297 bp
Locus ID:	25798
UniProt ID:	O95415
Cytogenetics:	7q21.3



Protein Families: Transmembrane

MW: 10.9 kDa

Gene Summary: Participates in tumor necrosis factor-alpha (TNF)-induced cell death (PubMed:14592447). May be a target of Wnt/beta-catenin signaling in the liver (PubMed:20538055).[UniProtKB/Swiss-Prot Function]