

Product datasheet for RC228057

C22orf39 (NM_001166242) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C22orf39 (NM_001166242) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C22orf39
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC228057 representing NM_001166242
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTGTCGGTGCTCATTAGTTCTTCTGTCAGTCGATCAGGAGTGCCTTTTCTTCATTTTTTCATTGGCT
 GGCGGACTGAGGGGCGGGCTTGGCGCGCCGTCGCCAGACATGGCGGACGGCAGCGGCTGGCAGCCGCC
 GCGCCCTGCGAGGCCTACCGCGCGAGTGGAGCTCTGCCGAGCGCCAGGCCTTCTACACCACTAC
 TACGTCCACGGCGAGCGCCGGCCTGCGAACAGTGGCAGCGCGACCTGGCCAGCTGCCGCGACTGGGAGG
 AGCGCCGGAACGCCGAGGCCAGGTGCTAAGTGTGGTGGCTCACACCTGTAATCTCAGCTACTTGGGAGG
 C

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228057 representing NM_001166242
 Red=Cloning site Green=Tags(s)

MCRCSLVLLSVDHEVPFSSFFIGWRTEGRAWRAGRPDMADGSGWQPPRPCEAYRAEWKLCRSARHFLHHY
 YVHGERPACEQWQRDLASCRDWEERRNAEQVLSVVAHTCNLSYLGG


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001166242

ORF Size: 351 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166242.1</u> , <u>NP_001159714.1</u>
RefSeq ORF:	243 bp
Locus ID:	128977
UniProt ID:	<u>Q6P5X5</u>
Cytogenetics:	22q11.21
MW:	13.4 kDa