

Product datasheet for **RC228064L4V**

zinc finger protein 138 (ZNF138) (NM_001160183) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	zinc finger protein 138 (ZNF138) (NM_001160183) Human Tagged ORF Clone Lentiviral Particle
Symbol:	zinc finger protein 138
Synonyms:	pHZ-32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001160183
ORF Size:	378 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228064).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001160183.1
RefSeq ORF:	381 bp
Locus ID:	7697
UniProt ID:	P52744
Cytogenetics:	7q11.21
Protein Families:	Transcription Factors
MW:	14.4 kDa


[View online »](#)

Gene Summary: May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]