

## Product datasheet for **RC228065L4V**

### ABCB5 (NM\_001163993) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ABCB5 (NM_001163993) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABCB5
Synonyms:	ABCB5alpha; ABCB5beta; EST422562
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163993
ORF Size:	378 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163993.1</a>
RefSeq Size:	1724 bp
RefSeq ORF:	381 bp
Locus ID:	340273
UniProt ID:	<a href="#">Q2M3G0</a>
Cytogenetics:	7p21.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ABC transporters


[View online »](#)

**MW:** 14.5 kDa

**Gene Summary:** ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules (Chen et al., 2005 [PubMed 15760339]).[supplied by OMIM, Mar 2008]