

## Product datasheet for **RC228075L3V**

### **HIST2H2BF (H2BC18) (NM\_001161334) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HIST2H2BF (H2BC18) (NM_001161334) Human Tagged ORF Clone Lentiviral Particle
Symbol:	H2BC18
Synonyms:	HIST2H2BF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001161334
ORF Size:	402 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228075).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161334.1</a> , <a href="#">NP_001154806.1</a>
RefSeq ORF:	405 bp
Locus ID:	440689
UniProt ID:	<a href="#">Q5QNW6</a>
Cytogenetics:	1q21.2
Protein Pathways:	Systemic lupus erythematosus
MW:	14.7 kDa



[View online »](#)

**Gene Summary:**

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-dependent histone that is a member of the histone H2B family and is found in a histone cluster on chromosome 1. [provided by RefSeq, Aug 2015]