

## Product datasheet for **RC228077L3V**

### ZNF268 (NM\_001165883) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF268 (NM_001165883) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF268
Synonyms:	HZF3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001165883
ORF Size:	405 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228077).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001165883.1</a>
RefSeq ORF:	408 bp
Locus ID:	10795
UniProt ID:	<a href="#">Q14587</a>
Cytogenetics:	12q24.33
Protein Families:	Transcription Factors
MW:	15.6 kDa



[View online »](#)

**Gene Summary:**

Isoform 1: Acts as a transcriptional repressor. Inhibits erythroid differentiation and tumor cell proliferation. Plays a role during ovarian cancer development and progression.  
[UniProtKB/Swiss-Prot Function]