

## Product datasheet for **RC228083**

### ZDHHC15 (NM\_001146257) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZDHHC15 (NM\_001146257) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ZDHHC15  
**Synonyms:** DHHC15; MRX91  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC228083 representing NM\_001146257  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

**ATGCGCGGAGGCTGGAAGATGGCTCTGTCTGGGGGCTGCGGTGCTGCCGCGGGTACTGTCTGGGTGC**  
**CAGTGCTCGTTATTGTCCTCGTCGTCTCTGGTCTACTATGCCTACGTCTTTGAACTCTGCCTGGTTAT**  
**TTACCTCATACTCTACCATGCCATCTTTGTGTTCTTTACCTGGACCTACTGGAAGTCTATCTTTACTC**  
**CCACAGCAGCCAAACCAGAAGTTCCACTTGTCTACACAGACAAGGAGCGCTATGAAAATGAAGAAAGAC**  
**CTGAGGTCAGAAAGCAGATGCTTGTGATATGGCCAAAAGCTACCGGTTTACACAAGAAGTGAAGTGG**  
**AGGTCAAGTTCATCAAAGGCAGCTAGAGAGGCAGCTCAGCAAGTATCTCAGAAAGGCTAAGTCATATATG**  
**TTCTCAAAC**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT**  
**ACAAGGATGACGACGATAAGGTTTAA**

**Protein Sequence:** >RC228083 representing NM\_001146257  
**Red=Cloning site Green=Tags(s)**  
  
MRRGWKMALSGGLRCCRVL SWVPVLVIVLVVLSYAYVFELCLVIYLILYHAIFVFFWTYWKSIFTL  
PQQPNQK FHL SYDKERYENEERPEVQKQMLVDMAKKLPVYTRTGSGGQFIQRQLERQLSKYLRKAKSYM  
FSN


**TRTRPLEQKLI SEEDLAANDILDYKDDDDKV**

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


**ACCN:** NM\_001146257

**ORF Size:** 429 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001146257.2</u>
<b>RefSeq ORF:</b>	432 bp
<b>Locus ID:</b>	158866
<b>UniProt ID:</b>	<u>Q96MV8</u>
<b>Cytogenetics:</b>	Xq13.3
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	17.1 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]