

## Product datasheet for **RC228088L4V**

### H2BFM (H2BW2) (NM\_001164416) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	H2BFM (H2BW2) (NM_001164416) Human Tagged ORF Clone Lentiviral Particle
Symbol:	H2BW2
Synonyms:	H2BFM; H2BM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001164416
ORF Size:	462 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228088).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164416.1</a> , <a href="#">NP_001157888.1</a>
RefSeq ORF:	447 bp
Locus ID:	286436
UniProt ID:	<a href="#">P0C1H6</a>
Cytogenetics:	Xq22.2
Protein Pathways:	Systemic lupus erythematosus
MW:	16.8 kDa



[View online »](#)

**Gene Summary:**

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene encodes a replication-independent histone that is a member of the H2B histone family. [provided by RefSeq, Nov 2015]