

Product datasheet for **RC228179A1V**

Human Endomucin (EMCN) (NM_001159694) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Endomucin (EMCN) (NM_001159694) AAV Particle
Tag:	Myc-DDK
Symbol:	Endomucin
Synonyms:	EMCN2; MUC14
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC228179 representing NM_001159694 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGCTGCTCAAGTGACCATTCCTTTTCTCTGCCAGTATTTGCAGCAGTAACAGCACAGGTGTTT
TAGAGGCAGCTAATAATTCACCTGTTGTTACTACAACAAAACCATCTATAACAACACCAACACAGAATC
ATTACAGAAAAATGTTGTCACACCAACAACCTGGAACAACCTCTAAAGGAACAATCACAATGAATTACTT
AAAATGTCTCTGATGTCAACAGCTACTTTTTTAACAAGTAAAGATGAAGGATTGAAAGCCACAACCACTG
ATGTCAGGAAGAATGACTCCATCATTTCAAACGTAACAGTAACAAGTGTTACACTTCCAAATGCTGTTTC
AACATTACAAAGTTCCAAACCAAGAGTAGTGTTCTACAACCAGATGCATCACCTTCTAAACTGGTACA
TTAACCTCAATACCAGTTACAATTCCAGAAAAACCTCACAGTCTCAAGTAATAGGCACTGAGGGTGGA
AAAATGCAAGCACTTCAGCAACCAGCCGGTCTTATTCCAGTATTATTTGCCGGTGGTTATTGCTTTGAT
TGTAATAACACTTTCAGTATTGTTCTGGTGGGTTTGTACCGAATGTGCTGGAAGGCAGATCCGGGCACA
CCAGAAAATGGAATGATCAACCTCAGTCTGATAAAGAGAGCGTGAAGCTTCTTACCGTTAAGACAATTT
CTCATGAGTCTGGTGAGCACTCTGCACAAGGAAAAACCAAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228179 representing NM_001159694 Red=Cloning site Green=Tags(s) MELLQVTILFLLPSICSSNSTGVLEAANNLSLVTTTKPSITTPNTESLQKNVVTPTTGTPKGTITNELL KMSLMSTATFLTSKDEGLKATTTDVRKNDISIISNVTVTSVTLPLNAVSTLQSSKPKSSVLQPDASPSKTGT LTSIPVTIPENTSQSQVIGTEGGKNASTSATRSYSSIILPVVIALIVITLSVFLVGLYRMCWKADPGT PENGNDQPQSDKESVKLLTVKTISHESGEHSAQGKTKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001159694
ORF Size:	744 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001159694.1</u>
RefSeq ORF:	747 bp
Locus ID:	51705
UniProt ID:	<u>Q9ULC0</u>
Cytogenetics:	4q24
MW:	25.9 kDa