

Product datasheet for **RC228196L4V**

CYB561A3 (NM_001161454) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CYB561A3 (NM_001161454) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CYB561A3
Synonyms:	CYBASC3; LCYTB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161454
ORF Size:	777 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001161454.1 , NP_001154926.1
RefSeq ORF:	780 bp
Locus ID:	220002
UniProt ID:	Q8NBI2
Cytogenetics:	11q12.2
Protein Families:	Transmembrane
MW:	29 kDa


[View online »](#)

Gene Summary:

Ferric-chelate reductase that reduces Fe(3+) to Fe(2+) before its transport from the endosome to the cytoplasm. Probably uses ascorbate as electron donor (By similarity).[UniProtKB/Swiss-Prot Function]