

Product datasheet for **RC228201L3V**

C17orf99 (NM_001163075) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C17orf99
Synonyms:	IL-40; UNQ464
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163075
ORF Size:	795 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228201).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001163075.1
RefSeq ORF:	798 bp
Locus ID:	100141515
UniProt ID:	Q6UX52
Cytogenetics:	17q25.3
MW:	28.9 kDa



Gene Summary:

Probable B cell-associated cytokine that plays a role in the regulation of humoral immune responses. Involved in lymphocyte B cell development and immunoglobulin/IgA production. [UniProtKB/Swiss-Prot Function]