

Product datasheet for **RC228243A1V**

Human FHL1 (NM_001159701) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FHL1 (NM_001159701) AAV Particle
Tag:	Myc-DDK
Symbol:	FHL1
Synonyms:	FCMSU; FHL-1; FHL1A; FHL1B; FLH1A; KYOT; RBMX1A; RBMX1B; SLIM; SLIM-1; SLIM1; SLIMMER; XMPMA
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC228243 representing NM_001159701 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACTTTCTATGTGGCCTCTCTGGCTCTGGAGCTAATTTGGATGCTGAGTAGCCCCGAGGTCCCTCCA
GCTACAAGGTGGGCACCATGGCGGAGAAGTTTGACTGCCACTACTGCAGGGATCCCTTGCAAGGGAAGAA
GTATGTGCAAAAGGATGGCCACCACTGCTGCCTGAAATGCTTTGACAAGTTCTGTGCCAACACCTGTGTG
GAATGCCGCAAGCCCATCGGTGCGGACTCCAAGGAGGTGCACTATAAGAACCGCTTCTGGCATGACACCT
GCTTCCGCTGTGCCAAGTGCTTCAACCCCTTGCCCAATGAGACCTTTGTGGCCAAGGACAACAAGATCCT
GTGCAACAAGTGCACTCGGGAGGACTCCCCAAGTGCAAGGGGTGCTTCAAGGCCATTGTGGCAGGA
GATCAAAACGTGGAGTACAAGGGGACCGTCTGGCACAAGACTGCTTCACTGTAGTAAGTCAAGCAAG
TCATCGGACTGGAAGCTTCTTCCCTAAAGGGGAGGACTTCTACTGCGTGACTTGCCATGAGACCAAGTT
TGCCAAGCATTGCGTGAAGTGCAACAAGGCCATCACATCTGGAGGAATCACTTACCAGGATCAGCCCTGG
CATGCCGATTGCTTTGTGTGTGTACCTGCTCTAAGAAGCTGGCTGGGCAGCGTTTACCAGCTGTGGAGG
ACCAAGTATTACTGCGTGGATTGCTACAAGAAGTTGTGGCCAAGAAGTGTGCTGGATGCAAGAACCCCAT
CACTGGGTTTGGTAAAGGCTCCAGTGTGGTGGCCTATGAAGGACAATCCTGGCACGACTACTGCTCCAC
TGCAAAAAATGCTCCGTAATCTGGCCAACAAGCGCTTTGTTTTCCACCAGGAGCAAGTATTGTCCTG
ACTGTGCCAAAAAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228243 representing NM_001159701 Red=Cloning site Green=Tags(s) MTFYVASLLELIWMLSSPAGPSSYKVGTMAEKFDCHYCRDPLQGKKYVQKDGHHCLKCFDKFCANTCV ECRKPIGADSKEVHYKNRFWHDTCFRCAKCLHPLANETFVAKDNKILCNKCTTREDSPKCKGCFKAIVAG DQNVVEYKGTVWHKDCFTCSNCKQVIGTGSFFPKGEDFYCVTCHETKFAKHCVKCNKAITSGGITYQDQPW HADCFVCVTCSKKLAGQRFTAVEDQYYCVDYKNFVAKKCAKCKNPITGFGKSSVVAIEGQSWHDYCFH CKKCSVNLANKRFVFHQEQVYCPDCAKKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001159701
ORF Size:	927 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001159701.1
RefSeq ORF:	930 bp
Locus ID:	2273
UniProt ID:	Q13642
Cytogenetics:	Xq26.3
MW:	34.8 kDa