

Product datasheet for **RC228262A1V**

Human RDH13 (NM_001145971) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RDH13 (NM_001145971) AAV Particle
Tag:	Myc-DDK
Symbol:	RDH13
Synonyms:	SDR7C3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC228262 representing NM_001145971 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCCGCTACCTGCTGCCGCTGTCCGCGCTGGGCACGGTAGCAGGCGCCGCCGTGCTGCTCAAGGACT
ATGTCACCGGTGGGGCTTGCCCCAGCAAGGCCACCATCCCTGGGAAGACGGTCATCGTGACGGGTGCCAA
CACAGGCATCGGAAGCAGACCGCCTTGGAAGTGGCCAGGAGAGGAGGCAACATCATCTGGCCTGCCGA
GACATGGAGAAGTGTGAGGCGGCAGCAAAGGACATCCGCGGGGAGACCCTCAATCACCATGTCAACGCCC
GGCACCTGGACTTGGCTTCCCTCAAGTCTATCCGAGAGTTTGCAGCAAAGATCATTGAAGAGGAGGAGCG
AGTGGACATTCTAATCAACAACGCGGGTGTGATGCGGTGCCCCCACTGGACCACCGAGGACGGCTTCGAG
ATGCAGTTTGGCGTTAACCACCTGGGTCACTTTCTCTTGACAACTTGCTGCTGGACAAGCTGAAAGCCT
CAGCCCCTTCGCGGATCATCAACCTCTCGTCCCTGGCCCATGTTGCTGGGCACATAGACTTTGACGACTT
GAACTGGCAGACGAGGAAGTATAACACCAAAGCCGCTACTGCCAGAGCAAGCTCGCCATCGTCTCTTC
ACCAAGGAGCTGAGCCGGCGGCTGCAAGGCTCTGGTGTGACTGTCAACGCCCTGCACCCCGGCGTGGCCA
GGACAGAGCTGGGCAGACACACGGGCATCCATGGCTCCACCTTCTCCAGCACCACACTCGGGCCCATCTT
CTGGCTGCTGGTCAAGAGCCCCGAGCTGGCCGCCAGCCAGCACATACCTGGCCGTGGCGGAGGAACTG
GCGGATGTTTCCGAAAGTACTTCGATGGACTCAAACAGAAGGCCCGGCCCGGAGGCTGAGGATGAGG
AGGTGGCCCGGAGGCTTTGGGCTGAAAGTGCCCGCCTGGTGGGCTTAGAGGCTCCCTCTGTGAGGGAGCA
GCCCCTCCCCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC228262 representing NM_001145971 Red=Cloning site Green=Tags(s) MSRYLLPLSALGTVAGAAVLLKDYVTGGACPSKATIPGKTVIVTGANTGIGKQTALARRGGNIILACR DMEKCEAAKDIRGETLNHHVNARHLDASLKSIREFAAKIIIEEEERVDILINNAGVMRCPHWTTEDGFE MQFGVNLGHFLLTNLLLDKLGASAPSRINLSSLAHVAGHIDFDDLNWQTRKYNTKAAYCQSKLAIVLF TKELSRRLQSGSVTVNALHPGVARTELGRHTGIHGSTFSSTTLGPFWLLVKSPELAAQPSTYLAVAEEL ADVSGKYFDGLKQKAPAEAEDEEVARRLWAESARLVGLEAPSVREQLPR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001145971
ORF Size:	993 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001145971.1</u>
RefSeq ORF:	996 bp
Locus ID:	112724
UniProt ID:	<u>Q8NBN7</u>
Cytogenetics:	19q13.42
MW:	35.8 kDa