

Product datasheet for **RC228268L3V**

MOV10L1 (NM_001164106) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	MOV10L1
Synonyms:	CHAMP; DJ402G11.8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164106
ORF Size:	1014 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228268).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001164106.1](#), [NP_001157578.1](#)

RefSeq ORF: 1017 bp

Locus ID: 54456

UniProt ID: [Q9BXT6](#)

Cytogenetics: 22q13.33

MW: 37.3 kDa



Gene Summary:

This gene is similar to a mouse gene that encodes a putative RNA helicase and shows testis-specific expression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]