

Product datasheet for **RC228269L4V**

NEK6 (NM_001166169) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NEK6
Synonyms:	SID6-1512
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166169
ORF Size:	1014 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228269).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001166169.1](#), [NP_001159641.1](#)

RefSeq ORF: 1017 bp

Locus ID: 10783

UniProt ID: [Q9HC98](#)

Cytogenetics: 9q33.3

Protein Families: Druggable Genome, Protein Kinase



MW: 38.3 kDa

Gene Summary: The protein encoded by this gene is a kinase required for progression through the metaphase portion of mitosis. Inhibition of the encoded protein can lead to apoptosis. This protein also can enhance tumorigenesis by suppressing tumor cell senescence. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]