

Product datasheet for **RC228282**

NAT1 (NM_001160176) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NAT1 (NM_001160176) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NAT1
Synonyms:	AAC1; MNAT; NAT-1; NATI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC228282 representing NM_001160176 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGTTATTACTCTTACACAAGGAGGCAGCCCTCGAGCCACAGGGTCCAGCTGTTGGCTATAATAGCC
TACCGGTCTCTGATGATCACCATGTTTCTGGAATCAAGCCAGGAAGAAGCAGCAATCTGTCTTCTGGAT
TAAACTGAAGATCAACCTACTTTCACTTAAGAAAGGGGATCATGGACATTGAAGCATATCTTGAA
AGAATTGGCTATAAGAAGTCTAGGAACAAATTGGACTTGGAAACATTAAGTACATTCTTCAACACCAGA
TCCGAGCTGTTCCCTTTGAGAACCTAACATCCATTGTGGGGATGCCATGGACTTAGGCTTAGAGGCCAT
TTTTGATCAAGTTGTGAGAAGAAATCGGGTGGATGGTGTCTCCAGGTCATCATCTTGTACTGGCT
CTGACCACTATTGGTTTTGAGACCACGATGTTGGGAGGGTATGTTTACAGCACTCCAGCCAAAAATACA
GCACTGGCATGATTCACCTTCTCCTGCAGGTGACCATTGATGGCAGGAACATATTGTCGATGCTGGGT
TGGACGCTCATACCAGATGTGGCAGCCTCTGGAGTTAATTTCTGGGAAGGATCAGCCTCAGGTGCCTTGT
GTCTTCCGTTTGACGGAAGAGAAATGGATTCTGGTATCTAGACCAATCAGAAGGGAACAGTACATTCCAA
ATGAAGAATTTCTTCTGATCTCTAGAAAGACAGCAATACCGAAAAATCTACTCCTTACTCTTAA
GCCTCGAACAATTGAAGATTTTGGTCTATGAATACATACCTGCAGACATCTCCATCATCTGTGTTTACT
AGTAAATCATTTTGTTCCTTGCAGACCCAGATGGGTTCACTGTTTGGTGGGTTCCACCTCACCATA
GGAGATTCAATTATAAGGACAATACAGATCTAATAGAGTTCAAGACTCTGAGTGAGGAAGAAATAGAAAA
AGTGCTGAAAAATATATTTAATATTTCTTGCAGAGAAAGCTTGTGCCAAACATGGTGATAGATTTTT
ACTATT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RC228282 representing NM_001160176
Red=Cloning site Green=Tags(s)

MLLLLLHKEAALEPQGPVAVGYNLSPVSDHHVSGIQARKKQQSVFWIKTEDQPTFNLLRKGIMDIEAYLE
 RIGYKKS RNKLDLETLTDILQHQIRAVPFENLNIHCGDAMDLEAIFDQVRRNRGGWCLQVNHLLYWA
 LTTIGFETMLGGVYVSTPAKKYSTGMIHLLQVTIDGRNYIVDAGFGRSYQMWPLELISGKDQPQVPC
 VFRLTEENGFWYLDQIRREQYIPNEEFHLSDDLLEDSKYRKIYSFTLKPRTIEDFESMNTYLQTSPPSVFT
 SKSFCSLQTPDGVHCLVGFVTLTHRRFNYKDNTDLIEFKTLSEEEIEKVLKNI FNISLQRKLVPKHGDRFF
 TI

TRTRPLEQKLISEEDLAANDILDYKDDDDKVV

Chromatograms: https://cdn.origene.com/chromatograms/mk8054_d05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN: NM_001160176

ORF Size: 1056 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160176.4](#)

RefSeq Size: 1939 bp

RefSeq ORF: 1059 bp

Locus ID: 9


Cytogenetics: 8p22

Protein Pathways: Caffeine metabolism, Drug metabolism - other enzymes, Metabolic pathways

MW: 40.8 kDa

Gene Summary: This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RC228282