

Product datasheet for RC228302

ER81 (ETV1) (NM_001163152) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ER81 (ETV1) (NM_001163152) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ER81
Synonyms: ER81
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC228302 representing NM_001163152
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTTCAAGATTTAAGTGCAAGTGCTTCTTTCCACCTGTTCACAACACAGAACGTTAGCTCAGGTAC
 CTGACAATGATGAGCAGTTTGTACCAGACTATCAGGCTGAAAGTTTGGCTTTTCATGGCCTGCCACTGAA
 AATCAAGAAAGAACCCACAGTCCATGTTCAAGAAATCAGCTCTGCCTGCAGTCAAGAACAGCCCTTTAAA
 TTCAGCTATGGAGAAAAGTGCCTGTACAATGTCAGATTTGCGCCGACGCTTTCTGAACCTGTAACCTCT
 TTCTCTTTTGGCAGCATGCCAAGGGAAGGACGTCCTATGTACCAACGCCAGATGTCTGAGCCAAACAT
 CCCCTTCCCACCACAAGGCTTTAAGCAGGAGTACCACGCCAGTGTATGAACACAACACCATGGTTGGC
 AGTGCGGCCAGCCAAAGCTTTCCCTCTCTGATGATTAACAGGAACCCAGAGATTTTGCATATGACT
 CAGAAGTGCTAGCTGCCACTCCATTTATATGAGGCAAGAAGGCTTCTGGCTCATCCCAGCAGAACAGA
 AGGCTGTATGTTTAAAAGGGCCCCAGGCAGTTTATGATGACACCTGTGTTGTCCAGAAAAATTCGAT
 GGAGACATCAAACAAGAGCCAGGAATGTATCGGAAGGACCCACATACCAACGGCGAGGATCACTTCAGC
 TCTGGCAGTTTTTGGTAGCTCTTCTGGATGACCCTTCAAATTCATTTTATTGCCTGGACTGGTCGAGG
 CATGGAATTTAACTGATTGAGCCTGAAGAGGTGGCCGACGTTGGGGCATTAGAAAAACAGGCCAGCT
 ATGAACTATGATAAACTTAGCCGTTCACTCCGCTATTACTATGAGAAAGGAATTATGCAAAAGGTGGCTG
 GAGAGAGATATGTCTACAAGTTTGTGTGTGATCCAGAAGCCCTTTCTCCATGGCCTTTCCAGATAATCA
 GCGTCCACTGCTGAAGACAGACATGGAACGTCACATCAACGAGGAGGACACAGTGCCTTTTCTCACTTT
 GATGAGAGCATGGCTACATGCCGGAAGGGGGCTGCTGCAACCCCAACCCTACAACGAAGGCTACGTGT
 AT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC228302 representing NM_001163152
Red=Cloning site Green=Tags(s)

MLQDLSASVFFPPCSQHRTLAQVPDNDQFVPDYQAESLAFHGLPLKIKKEPHSPCSEISSACSQEOPFK
 FSYGEKCLYNVRFRRQLSEPCNSFPPLPTMPREGRPMYQRMSEPNIFFPPQGFQYHDPVYEHTMVG
 SAASQSFPPPLMIKQEPDFAYDSEVPVSCHSIYMRQEGFLAHPSTRTEGCMFEKGPQRFYDDTCVPEKFD
 GDIKQEPGMYREGPTYQRRGSLQLWQFLVALLDDPSNSHFIAWTGRGMEFKLIEPEEVARRWGIQKNRPA
 MNYDKLSRSLRYYYEKIMQKVAGERVYVYKFCVDPALF SMAFPDNRPLLKTDMERHINEEDTVPLSHF
 DESMAYMPEGCCNPHPYNEGYVY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8055_e11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001163152

ORF Size: 1122 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001163152.1, NP_001156624.1</u>
RefSeq ORF:	1125 bp
Locus ID:	2115
UniProt ID:	<u>P50549</u>
Cytogenetics:	7p21.2
Protein Families:	ES Cell Differentiation/IPS, Transcription Factors
MW:	43.3 kDa
Gene Summary:	<p>This gene encodes a member of the ETS (E twenty-six) family of transcription factors. The ETS proteins regulate many target genes that modulate biological processes like cell growth, angiogenesis, migration, proliferation and differentiation. All ETS proteins contain an ETS DNA-binding domain that binds to DNA sequences containing the consensus 5'-CGGA[AT]-3'. The protein encoded by this gene contains a conserved short acidic transactivation domain (TAD) in the N-terminal region, in addition to the ETS DNA-binding domain in the C-terminal region. This gene is involved in chromosomal translocations, which result in multiple fusion proteins including EWS-ETV1 in Ewing sarcoma and at least 10 ETV1 partners (see PMID: 19657377, Table 1) in prostate cancer. In addition to chromosomal rearrangement, this gene is overexpressed in prostate cancer, melanoma and gastrointestinal stromal tumor. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2016]</p>

Product images:



Circular map for RC228302



Western blot validation of overexpression lysate (Cat# [LY431330]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC228302 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).