

Product datasheet for **RC228315L3V**

ZNF414 (NM_001146175) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZNF414
Synonyms:	ZFP414
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146175
ORF Size:	1170 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC228315).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001146175.1
RefSeq ORF:	1173 bp
Locus ID:	84330
UniProt ID:	Q96IQ9
Cytogenetics:	19p13.2
MW:	40.8 kDa



Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]