

Product datasheet for **RC228327**

ZNF735 (NM_001159524) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZNF735
Synonyms:	ZNF735P
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RC228327 ORF sequence, **codon optimized**.
Due to the complexity of NM_001159524, the ORF clone is codon optimized for mammalian Expression.
The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCAAAGAGGCCTGGACCCCGGGGAGCAGAGAAATGGCCTGCTGACGTTCCGGGACATCGCAATCG
AATTCAAGTTTGGCGGAGTGGCAATGCCTGGATCATGCCAGCAGAATCTTTATAGAGATGTGATGTTGGA
GAATATAGGAATCTCTTAGCCTGGGTATGACTGTGTCAAACCTGATCTTATTGCTTGCTCGAACAG
AATAAGGAGCCCCAAAACATCAAAAGAAACGAGATGGCTGCTAAACATCCTGTGACCTGTAGCCACTTCA
ACCAAGATCTGCAGCCAGAGCAGTCCATAAAGGACAGTTTGCAGAAAGTCATACCCCGCACATACGGGAA
GTGTGGGCACGAGAACCTTCAGTTGAAGAAGTGTGTAACGGGTGGACGAATGCGAGGTGCATAAGGGG
GGGTACAACGACCTGAATCAATGCCTCTCTAATACGCAGAACAAAATTTTTCAGACCCACAAATGCGTGA
AGGTATTCTCTAAGTTTAGTAATTCTAACCGACACAACGCCAGATACACTGGGAAGAAACATCTCAAATG
TAAGAAATACGGCAAGTCTTTTGTATGTTCTCCACCTGAACCGACACAGATCATTATACCAAGGAG
AAAAGTTATAAATGCGAGGAGTGGCGGAAGTCTTTAATCACTCTCTTCCGGGACCAACCCATAAGCGAA
TTCTGACGGGAGAAAAACCATACCGCTGTGAGGAATGCGGCAAGCGTTCGCTGGCCATCAAATCTCAC
CAGGCACAAAAGAAATTCACACAGGGGAGAACCCCTACGCTTGTGAGGAGTGGCGCAAGCCTTTAGACGC
TCCTCCACCTGACAAATCATAGCGGATTCTACAGGAGAGCGCCGTATAAATGCGAAGAGTGTGGGA
AGGCGTTCTCCGTTTCTCCGCACTGATCTACCATAAAGAATTACACAGGAGAGAGCCTTACACCTG
TGAGGAATGCGGGAAGGCCTTAACTGCAGCTCAACCTGAAGACGCATAAGATAATACATACCGGTGAG
AAACCTTATACCTGTGAAGAGTGGGTGCAACGTTCAATTGTAGTAGTACTGTAAAGGCACACAAGCGCA
TCCACACGGGCGAAAAAGCCCTACAAATGCGAAGAATGCGACAAAGCATTTAAGTGGCATAGCAGTCTGGC
CAAGCATAAAATTATCCACACCGGGGAAAAACCATATAAGTGCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228327 representing NM_001159524
Red=Cloning site Green=Tags(s)

MAKRPGPPGSREMGLLTFRDIAIEFSLAEWQCLDHAQQNLYRDVMLENYRNLFSLGMTVSKPDLIACLEQ
NKEPQNIKRNEMAAKHPVTCSHFNQDLQPEQSIKDSLQKVIPTRYGKCGHENLQLKKCKRVDECEVHK
GYNDLNQCLSNQNKIFQTHKCVKVFSEKFSNSNRHNARYTGKKHLKCKKYGKSFCMFSHLNQHQIHTKE
KSYKCEECGKSFNHSSSGTTHKRILTGEKPYRCEECGKAFRWPSNLTRHKRIHTGEKPYACEECGQAFRR
SSTLTNHNKRIHTGERPKCEECGKAFSVSSALIYHKRIHTGEKPYTCEECGKAFNCSSTLTKHKIHTGE
KPYTCEECGRFNCSSSTVKAHKRIHTGEKPYKCEECDAFKWHSSSLAKHKIHTGEKPYKCK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001159524

ORF Size: 1236 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001159524.1](#), [NP_001152996.1](#)

RefSeq Size: 1239 bp

RefSeq ORF: 1239 bp

Locus ID: 730291


UniProt ID: [P0CB33](#)

Cytogenetics: 7q11.21

MW: 47.6 kDa

Gene Summary: This gene encodes a kruppel-associated box-containing zinc finger protein (KRAB-ZFP). The encoded protein contains an N-terminal kruppel-associated box (KRAB) domain and nine C-terminal C2H2-type zinc finger domains. The KRAB-ZFPs represent the largest family of mammalian transcriptional repressors, which function through the recruitment of the nuclear co-factor KRAB-Associated Protein 1 (KAP1), to engage histone modifiers and induce heterochromatin formation. [provided by RefSeq, Jul 2017]

Product images:



Circular map for RC228327