

## Product datasheet for **RC228356L3V**

### SLC14A1 (NM\_001146037) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC14A1 (NM_001146037) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC14A1
Synonyms:	HsT1341; HUT11; JK; Jk(b); RACH1; RACH2; UT-B1; UT1; UTE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146037
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228356).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001146037.1</a> , <a href="#">NP_001139509.1</a>
RefSeq ORF:	1338 bp
Locus ID:	6563
UniProt ID:	<a href="#">Q13336</a>
Cytogenetics:	18q12.3
Protein Families:	Transmembrane
MW:	48.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a membrane transporter that mediates urea transport in erythrocytes. This gene forms the basis for the Kidd blood group system. [provided by RefSeq, Mar 2009]