

Product datasheet for **RC228379**

FAM90A9 (NM_001164448) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM90A9
Synonyms:	FAM90A9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC228379 representing NM_001164448
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGATGGCAGTCGGGACCCCAATCTTGGCCAAGAGACTGGTGAGAGCCAGACCCTCCAGAAGCAGC
 GGAGGGCCCCAGTTGGGCCAAGGGCTCCCCCGCCGATGAAGAAGATCCCAGGCTCAAGTGCAAAAACCTG
 CGGGGCTTTGGCCACACGGCCAGAAGTACCAGGTGCCCCATGAAGTGCTGGAAGGCAGCCCTGGTTCCA
 GCGACCTTGGGGAAAAAGGAAGGGAAGGAAACCTGAAACCATGGAAGCCCCGGGTTGAAGCCAACCCGG
 GGCCCTTGAACAAGGATAAGGGAGAGAAGGAAGAGAGACCAAGGCAACAAGACCCGAGAGGAAGGCTCT
 CCTCCACATGTTTTCCGGGAAACCTCCAGAGAAGCCGCTGCCGAATGGAAAAGGATCCACGGAGTCTTCT
 GATCATCTGAGGTTGCAAGCGGGCAATGCCGGTCCACACAACCAAGTAAGAGGCCGCGCTGGACCCTG
 TCCTCGCTGATCGCTCAGCTGCCGAAATGTCTGGCAGGGGCTCCGTCTTGGCTTCACTGTCTCCCTCAG
 AAAAGCCAGCCTGAGCTCCTCTCAAGTCTTGGACCAAGGAAAGACAGACAGGGGCTGCGGCCGACATG
 CCTCAGCCTGCAGTCAGGCACAGGGCCGCGAGCCTCTCCTCGTGGTGAAGCCGACACACAGCAGCCCCG
 AGGGTGGCTGCCGAGAAGTTCCCCAGGCTGCCTCCAAAACCCACGGCCTGCTCCAGGCCGCCAGACCCCA
 GGCACAAGACAAACGTCTCGGGTACCTCGCAGCCCTGCCGCCAGCCGCCACACAGCTTGGGCCTA
 GGCTCAATCTCAGCTTCGGGCCAGGAGCCAAGAGACTGCCAGGCTCCGATTCAGGCTTGCCTGAAC
 TCCCCAAGAAACCGAGACTGGGTCCCTTCCAGATCCCCGAAAGCGCCATCCAGGGAGGTGAGCTGGGGGC
 CCCGGAGAATCTCAACCTCCGCCAGCCGCAACCGAACTTGACCAAGTACGTGCCCCAGATGGGCAGG
 AGGACACCGGCCAGGTGCCAGCGTCGACCGGCAGCCTCCGCACAGCACACCTTGCCTGCCTACTGCC
 AGGCCTGCACCATGTCCATCACTCAGCGGCCAGCCATGATGGGGCCAGCCTCTCAGAGTGCTCTTCCG
 GAGACTGGAACCGAGCCTGGAGCTCCAGCCTCCTGGCGGCCCTCATTTCACTCTCTGAGAAGCCG
 GGAACCTTCTCGCTCAGAGCCTCATGTGTGAGAGAAGTCTGAGGCTCCCTGTGTGCTGCCACCGA
 GCGTCTCTATGAGGACCTCAGGTTCTCTCTCCTCAGAGGACAGCGATTCTGACCTGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228379 representing NM_001164448
 Red=Cloning site Green=Tags(s)

MMARRDPKSWAKRLVRAQTLQKQRRAPVGPRAAPPDEEDPRLKCKNCGAFGHTARSTRCPMKCWKAALVP
 ATLKGKEGKENLKPWKPRVEANPGPLNKDKEKEERPRQDPQRKALLHMFSGKPPEKPLPNGKSTESS
 DHLRVASGPMVHTTSKRPRVDPVLADRSAAEMSGRGSVLASLSPLRKASLSSSSSLGPKERQTGAADM
 PQPAVRHQGREPLL VVKPTHSSPEGGCREVPQAASKTHGLLQAARPAQDKRPVTSQPCPPAATHSLGL
 GSNLSFGPGAQRPAQAPIQAACLNFPKPKRLGPFQIPESAIQGGELGAPENLQPPPAATELGPSTSPQMGR
 RTPAQVPSVDRQPPHSTPCLPTAQAQCTMSHSAASHDGAQPLRVLFRRLENGRWSSSLLAAPSFSHSEPKP
 GTFLAQSPHVSEKSEAPCVRVPPSVLYEDLQVSSSSSDSDSLE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8038_g04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001164448

ORF Size: 1392 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001164448.1](#), [NP_001157920.1](#)

RefSeq ORF: 1394 bp


Locus ID: 441327

Cytogenetics: 8p23.1


MW: 49.6 kDa

Gene Summary: FAM90A9 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, Oct 2009]


Product images:



Circular map for RC228379



Western blot validation of overexpression lysate (Cat# [LY431407]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC228379 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified FAM90A9 protein (Cat# [TP328379]). The protein was produced from HEK293T cells transfected with FAM90A9 cDNA clone (Cat# RC228379) using MegaTran 2.0 (Cat# [TT210002]).