

Product datasheet for **RC228385**

FAM90A5 (NM_001164455) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM90A5
Synonyms:	FAM90A5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC228385 representing NM_001164455
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATGGCAGTCGGGACCCACATCTTGGCCAAGAGACTGGTGAGAGCCAGACCCTCCAGAAGCAGC
 GGAGGGCCCCAGTTGGGCCAAGGGCTCCCCCGCCGATGAAGAAGATCCCAGGCTCAAGTGCAAAAACCTG
 CGGGGCTTTGGCCACACGGCCAGAAGTACCAGGTGCCCATGAAGTGCTGGAAGGCAGCCCTGTTTCCA
 GCGACCTTGGGGAAAAAGGAAGGGAAGGAAACCTGAAACCATGGAAGCCCCGGGTGAAGCCAACCCGG
 GGCCCTTGAACAAGGATAAGGGAGAGAAGGAAGAGAGACCAAGGCAACAAGACCCGAGAGGAAGGCTCT
 CCTCCACATGTTTTCCGGGAAACCTCCAGAGAAGCCGCTGCCGAATGGAAAAGGATCCACGGAACCTTCT
 GATTATCTGAGGTTGCAAGCGGGCAATGCCGGTCCACACAACCAAGTAAGAGGCCGCGCTGGACCCTG
 TCCTCGCTGATGGCTCAGCTACCGAAATGTCTGACAGGGGCTCCGTCTTGGCTTCACTGTCTCCCTCAG
 AAAAGCCAGCCTGAGCTCCTCTCAAGTCTTGGACCAAGGAAAGACAGACAGGGGCTGCGGCCGACATG
 CCTCAGCCTGCAGTCAGGCACAGGGCCGCGAGCCTCTCCTCGTGGTGAAGCCGACACACAGCCGCCCG
 AGGGTGGCTGCCGAGAAGTTCCCCAGGCTGCCTCCAAAACCCACGGCCTGCTCCAGGCCGCCAGACCCCA
 GGCACAAGACAAACGTCTCGGGTACCTCACAGCCCTGCCCGCCAGCCGCCACACAGCTTGGGCCTA
 GGCTCAATCTCAGCTTCGGGCCAGGAGCCAAGAGACTGCCAGGCTCCGATTCAGGCTTGCTGAACT
 TCCCCAAGAAACCGAGACTGGGTCCCTTCCAGATCCCCGAAAGCGCCATCCAGGGAGGTGAGCTGGGGGC
 CCCGGAGAATCTCAACCTCCGCCAGCCGCAACCGAACTTGACCAAGTACGTGCCCCAGATGGGCAGG
 AGGACACCGGCCAGGTGCCAGCGTCGACCGGCAGCCTCCGCACAGCAGACCTTGCTGCCACTGCCC
 AGGCTGCACCATGTCCATCAGCCAGCGGCCAGCCATGATGGGGCCAGCCTCTCAGAGTGCTCTTCCG
 GAGACTGGAACCGAGCCTGGAGCTCCAGCCTCTGGCGGCCCTCATTTCACTCTCTGAGAAGCCG
 GGAGCCTTCTCGCTCAGAGCCTCATGTGTGAGAGAAGTCTGAGGCTCCCTGTGTGTCGCCACCGA
 GCGTCTCTATGAGGACCTCAGGTTCTCTCTCTCAGAGGACAGCGATTCTGACCTGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228385 representing NM_001164455
 Red=Cloning site Green=Tags(s)


MMARRDPTSWAKRLVRAQTLQKQRRAPVGPRAAPPDEEDPRLKCKNCGAFGHTARSTRCPMKCWAALVP
 ATLKGKEGKENLKPWKPRVEANPGPLNKDKEKEERPRQDPQRKALLHMFSGKPPEKPLPNGKSTEPS
 DYLRVASGPMVHTTSKRPRVDPVLADGSATEMSDRGSVLASLSPLRKASLSSSSSLGPKERQTGAADM
 PQPAVRHQGREPLL VVKPTHSRPEGGCREVPQAASKTHGLLQAARPAQDKRPVTSQPCPPAATHSLGL
 GSNLSFGPGAKRPAQAPIQAQLNFPKPKRLGPFQIPESAIQGGELGAPENLQPPPAATELGPSTSPQMGR
 RTPAQVPSVDRQPPHSRPLPTAQAQCTMSHHPAASHDGAQPLRVLFRRLENGRWSSLLAAPSFSHSEPK
 GAFLAQSPHVSEKSEAPCVRVPPSVLYEDLQVSSSSSDSDSLE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8038_d12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001164455

ORF Size: 1392 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22µm filter is required.

RefSeq: [NM_001164455.1](#), [NP_001157927.1](#)

RefSeq ORF: 1394 bp


Locus ID: 441315

Cytogenetics: 8p23.1

MW: 49.7 kDa


Gene Summary: FAM90A5 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, May 2010]

Product images:



Circular map for RC228385

Western validation with an anti-DDK antibody
* L: Control HEK293 lysate R: Over-expression lysate



Coomassie blue staining of purified FAM90A5 protein (Cat# [TP328385]). The protein was produced from HEK293T cells transfected with FAM90A5 cDNA clone (Cat# RC228385) using MegaTran 2.0 (Cat# [TT210002]).