

Product datasheet for **RC228406**

MTF2 (NM_001164391) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MTF2 (NM_001164391) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MTF2
Synonyms:	dj976O13.2; M96; PCL2; TDRD19A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC228406 representing NM_001164391
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGTCTGTACAATATGTCAAGAAGAGTATTCAGAAGCTCCCAATGAAATGGTTATATGTGACAAGTGTG
 GCCAAGGATATCATCAGTTGTGTACACACCTCATATTGATTCCAGTGTGATTGATTGAGATGAAAAATG
 GCTCTGTGCGCAGTGTGTTTTGCAACAACAACAAGAGGGGTGGTGCACCTTAAGAAAGGACCAAAATGCC
 AAAGCATTGCAAGTCATGAAGCAGACATTACCCTATAGTGTGGCAGACCTTGAATGGGATGCAGGTCATA
 AAACCAATGTCCAGCAGTGTACTGCTATTGTGGAGGCCCTGGAGACTGGTATTTGAAGATGCTACAGTG
 CTGCAATGTAAGCAGTGGTTTCATGAGGCTTGTGTGCAATGCCTTCAAAAGCCAATGCTATTTGGAGAC
 AGATTTTATACGTTTATATGCTCTGTCTGCAGTCTGGACCAGAATACCTCAAACGTCTACCATTACAGT
 GGGTAGATATAGCACACCTATGCCTTTACAACCTAAGTGTATTTCATAAGAAGAAATACTTTGATTCTGA
 ACTTGAGCTTATGACATACATTAATGAAAACCTGGGATAGATTGCACCCTGGAGAGCTGGCAGACACACCA
 AAATCTGAAAGATATGAGCATGTTCTGGAGGCATTAATGATTACAAGACCATGTTTATGCTGGGAAAG
 AAATAAAGAAGAAGAAGCATTGTTTGGGTTGCGAATTCGTGTTCTCTCTGTGCCACCAAAATGTGGCTTT
 CAAAGCAGAGAAAGAACCTGAAGGAACATCTCATGAATTTAAATTAAGGCAGAAAGGCATCCAAACCT
 ATATCTGATTCAGGGAAGTAAGCAATGGCATAGAAAAAAAGGAAAGAAAAATCTGTAGGTCGTCAC
 CTGGCCCATATACAAGAAAAATGATTCAAAAACTGCTGAGCCACTTTTGGAAGGAATCAATTTGAGA
 GAATCTACTTTGGATTTACCTTGTTCTATAGGGAGAAGTGGGAACTGCACATTCATCAATACCTCA
 GATGTGGATTTACGGGTGCTTCCAGTGCAAAAGAAACTACCTCGTCTAGCATTTCAGGCATTATGGAT
 TATCTGACTCCAGAAAAAGACGCGTACAGGAAGATCTTGGCCTGCTGCAATACCAATTCGCGGGAAG
 AAGAGGTCGTCTTCCAAGAAGAGCACTCCAGACTCAGAAGTGTAAAAGATGATGAAGGCAAA
 GAAGATTACAGTTTGTGAAGTCAACACAGATTCTGAATAACTTAGCAGATCAGGAGTTACAACCTCA
 ATCATCTAAAGAACTCCATTACCAGTTATTTGGTGTGTCAGGTAGAATAGCATGTGGCGAAAAATACCG
 AGTTTTGGCACGTCGGGTGACACTTGATGGAAAGTGCAGTATCTTGTGGAATGGGAAGGAGCAACTGCA
 TCC

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC228406 representing NM_001164391
 Red=Cloning site Green=Tags(s)

MVCTICQEEYSEAPNEMVICDKCGQGYHQLCHTPHIDSSVIDSDEKWLCRQCVFATTTKRGGALKKGPNA
 KALQVMKQTL PYSVADLEWDAGHKTNVQQCYCYCGGPGDWYLMKLQCKCKQWFHEACVQLQKPM LFGD
 RFYTFICSVCSGPEYLKRLPLQWVDIAHLCLYNLSVIHKKKYFDSELELMTYINENWDRLHPGELADTP
 KSERYEHVLEALNDYKTMFMSGKEIKKKKHLFGLRIRVPPVPPNVAFKAEKEPEGTSHEFKIKGRKASKP
 ISDSREVSNGIEKKKKKSVGRPPGYPYTRKMIQKTAEP LLDKESISENPTLDLPCSIGRTEGTAHSSNTS
 DVDFTGASSAKETTSSSISRHYGLSDSRKRTRTGRSWPAAIPHLRRRRGRLPRRALQTQNSEIVKDDEGK
 EDYQFDELNTEILNNLADQELQLNHLKNSITSYFGAAGRIACGEKYRVLARRVTLDGKVQYLVEWEGATA
 S

SGP**TRRRLEQKLI**SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8063_b08.zip

Restriction Sites:

SgfI-RsrII

Cloning Scheme:


ACCN: NM_001164391

ORF Size: 1473 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164391.1](#), [NP_001157863.1](#)

RefSeq ORF: 1476 bp

Locus ID: 22823

UniProt ID: [Q9Y483](#)


Cytogenetics: 1p22.1

Protein Families: Druggable Genome, Transcription Factors

MW: 55.5 kDa

Gene Summary: Polycomb group (PcG) that specifically binds histone H3 trimethylated at 'Lys-36' (H3K36me3) and recruits the PRC2 complex. Acts by binding to H3K36me3, a mark for transcriptional activation, and recruiting the PRC2 complex, leading to enhance PRC2 H3K27me3 methylation activity. Regulates the transcriptional networks during embryonic stem cell self-renewal and differentiation. Promotes recruitment of the PRC2 complex to the inactive X chromosome in differentiating XX ES cells and PRC2 recruitment to target genes in undifferentiated ES cells. Required to repress Hox genes by enhancing H3K27me3 methylation of the PRC2 complex. In some conditions may act as an inhibitor of PRC2 activity: able to activate the CDKN2A gene and promote cellular senescence by suppressing the catalytic activity of the PRC2 complex locally. Binds to the metal-regulating-element (MRE) of MT1A gene promoter (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC228406