

## Product datasheet for **RC228484L4V**

### Synaptotagmin 14 (SYT14) (NM\_001146264) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Synaptotagmin 14 (SYT14) (NM_001146264) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Synaptotagmin 14
Synonyms:	SCAR11; sytXIV
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001146264
ORF Size:	1800 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228484).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001146264.1</a>
RefSeq ORF:	1803 bp
Locus ID:	255928
Cytogenetics:	1q32.2
Protein Families:	Transmembrane
MW:	67 kDa


[View online »](#)

**Gene Summary:**

This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that mediate membrane trafficking in synaptic transmission. The encoded protein is a calcium-independent synaptotagmin. Mutations in this gene are a cause of autosomal recessive spinocerebellar ataxia-11 (SCAR11), and a t(1;3) translocation of this gene has been associated with neurodevelopmental abnormalities. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Dec 2011]