

## Product datasheet for **RC228506L4V**

### Neurofascin (NFASC) (NM\_001005389) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Neurofascin (NFASC) (NM_001005389) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Neurofascin  |
| Synonyms:                 | NEDCPMD; NF; NRCAML  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001005389   |
| ORF Size:                 | 1857 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC228506).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001005389.1</a>   |
| RefSeq ORF:               | 1860 bp  |
| Locus ID:                 | 23114  |
| UniProt ID:               | <a href="#">O94856</a>   |
| Cytogenetics:             | 1q32.1   |
| Protein Families:         | Transmembrane  |
| Protein Pathways:         | Cell adhesion molecules (CAMs)   |
| MW:                       | 69.57 kDa  |



[View online »](#)

**Gene Summary:**

This gene encodes an L1 family immunoglobulin cell adhesion molecule with multiple IGcam and fibronectin domains. The protein functions in neurite outgrowth, neurite fasciculation, and organization of the axon initial segment (AIS) and nodes of Ranvier on axons during early development. Both the AIS and nodes of Ranvier contain high densities of voltage-gated Na<sup>+</sup> (Nav) channels which are clustered by interactions with cytoskeletal and scaffolding proteins including this protein, gliomedin, ankyrin 3 (ankyrin-G), and betaIV spectrin. This protein links the AIS extracellular matrix to the intracellular cytoskeleton. This gene undergoes extensive alternative splicing, and the full-length nature of some variants has not been determined. [provided by RefSeq, May 2009]