

Product datasheet for **RC228513**

ARHGAP25 (NM_001166276) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP25 (NM_001166276) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ARHGAP25
Synonyms:	HEL-S-308; KAIA0053
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC228513 representing NM_001166276
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCCCTCGGTCAAGGAGTGTGATGACTGGCGAGC
 AGATGGCTGCCTCCATCCATCGTCCACCCCAACCCGCTGGAGAGGCCATCAAGATGGGCTGGCTGAA
 GAAGCAGAGGTCCATCGTGAAGAACTGGCAGCAGAGTACTTTGTGCTGAGGGCGCAGCAGCTCTACTAC
 TACAAGGATGAAGAGGACACGAAGCCCCAGGGCTGCATGTATCTACCAGGATGTACAATCAAGGAGATCG
 CCACAAACCCAGAAGAAGCTGGGAAGTTTGTCTTTGAAATCATTCCAGCCTCATGGGACCAGAATCGCAT
 GGGACAGGACTCCTATGTCTCATGGCCAGCTCTCAGGCGGAGATGGAGGAGTGGTTAAATTCCTCAGG
 AGAGTTGTGGCACACCCTGTGGAGCAGTGTGGCCAGCGCTTGGATGAGACTGTGGCCTATGAACAGA
 AATTCGGCCCCCATCTGGTGCCCATCTGGTGGAGAAATGTGCAGAGTTCATCTGGAGCACGGCCGGAA
 TGAAGAGGGCATCTTCCGTCTGCCTGGCAGGACAACCTGGTGAAGCAGCTGAGAGACGCTTTTGATGCT
 GGGGAGCGGCCCTCCTTTGACAGAGACACAGATGTGCACACTGTGGCTTCCCTGTAAAGCTCTACCTCC
 GAGACCTCCCAGAGCCCCTGGTCCCTGGAGCCAGTACGAAGGGTTCTGCTCTGTGGGCAGCTCAGGAA
 TGCGGATGAGGCAAAGGCTCAGCAGGAGTTGATGAAGCAGCTCTCCATCCTTCTCGTGACAATAAGT
 CTCCTGAGCTACATCTGCAGGTTCTACATGAAATACAGCTGAACTGTGCTGTTAACAAGATGAGTGTGG
 ACAACCTGGTACTGTGATTGGTGTGAATCTCATCAGGTCGAAGGTCGAAGACCCCTGCCGTGATCATGAG
 AGGGACTCCTCAGATCCAAGAGTGTGACTATGATGATCAGAGACCATGAAGTCTCTTCCCAAGTCC
 AAGGATATACCCCTGTCACCCCTGCCAGAAAAATGACCCCAAGAAAGCTCCAGTGGCCCGAAGCTCTG
 TAGGCTGGGATGCCACTGAAGACCTCCGAATTTCTAGGACAGACAGCTTCAGTAGCATGACAAGCGACTC
 TGATAACAACAGCCCCACCGGACAGCAGCCGAGCGATGCGTTTCCGGAGGACAGCAGCAAAGTACCCAGG
 GAAAAGCCAGGAGACTGGAAAATGCAATCTCGTAAAAGGACTCAAACACTCCCTAACCGGAAATGTTTCT
 TGACATCAGCTTTTCAGGGTGCCAACAGCAGCAAAATGGAGATCTTAAAAATGAATTCGGTCGCCTTC
 CTCAGAGGCTAAGGCAGGGGAAGGGCACAGGAGAACGATGTCTCAAGACTTGGCCAACTTTCTGACTCC
 CAACGGACTTCCACCTACGATAACGTCCCTTCCCTGCCAGGGTCCCCTGGGGAGGAAGCCAGTGCCTCT
 CTTCCCAAGCCTGTGACTCCAAGGGAGATACTCTTCCAGTCCAACTCTGAACTGGGCCTGGAAAAAA
 GAACTCTGGAGAAGAGGAAATGATTCTTTGCAGAGGATGGTCCAAGAGCTACGAAAGGAAATAGAAACA
 CAGAAGCAAATGTATGAGGAACAGATTAACAACTTGGAGAAGGAAATATGACGTTTGGGCTAAAGTGG
 TGAGGCTCAATGAAGAACTGGAGAAGGAAAAGAAGTCTGCAGCCCTAGAGATCAGCCTCCGCAACAT
 GGAGCGCTCCCGGAGGATGTTGAGAAGAGGAACAAGGCCTTGGAGAAGAAGTCAAGGAATTTGTCAAA
 TCCATGAAGGAACCAAGACCGAGGCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228513 representing NM_001166276
 Red=Cloning site Green=Tags(s)


```
MSLGQSACLFLSIARSRVMTGEQMAAFHPSSTPNPLERPIKMGWLKKQRSIVKNWQQRVFLRAQQLYY
YKDEEDTKPQGCMYLPGCTIKEIATNPEEAGKFVFEIIPASWDQNRMGQDSYVLMASSQAEMEEWVKFLR
RVAGTPCGAVFGQRLDETVAYEQKFGPHLVPILVEKCAEFIEHGRNEEGIFRLPGQDNLVKQLRDAFDA
GERPSFDRDRTDVHTVASLLKLYLRDLPEPVPWSQYEGFLLCGQLTNADEAKAQQLMKQLSILPRDNYS
LLSYICRFLHEIQLNCAVNKMSVDNLATVIGVNLIRSKVEDPAVIMRGTPQIQRVMTMMIRDHEVLFPKS
KDIPLSPPAQKNDPKKAPVARSSVGDATEDLRI SRTDSFSSMTSDSDTT SPTGQQPSDAFPEDSSKVP
EKPGDWKMQSRKRTQTLPNRKCF LSAFQGANSSKMEIFKNEFWSPSSEAKAGEGHRRTMSQDLRQLSDS
QRTSTYDNVPSLPGSPGEEASALSSQACDSKGD L ASPNSETGPGKNSGEE EIDSLQRMVQELRKEIET
QKQMYEEQIKNLEKENYDVWAKVRLNEELEKEKKKSAALEISLRNMERSREDVEKRNKALEEEVKEFKV
SMKEPKTEA
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8071_a06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001166276

ORF Size: 1917 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166276.2](#)

RefSeq ORF: 1920 bp

Locus ID: 9938


UniProt ID: [P42331](#)

Cytogenetics: 2p13.3

MW: 72.3 kDa

Gene Summary: ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC228513