

Product datasheet for **RC228516**

ALS (IGFALS) (NM_001146006) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALS (IGFALS) (NM_001146006) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ALS
Synonyms:	ACLSL; ALS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC228516 representing NM_001146006
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCCTGAGGAAAGCGGGTGACCTCGAGCCTCAGTTACACCTGAGCGCCGCTCCGGCTGTGCTGGT
 ATCAGGCTCACTCGGGCAGAGCCCTCCTCGGGCCGCTCCCAAGCCTCTCCCCCGCAGGAGCCCTGGC
 CCTGGCGCTGCTGCTGCTCCTGGGTGGCACTGGGCCCCCGCAGCCTGGAGGGAGCAGACCCCGGAACG
 CCGGGGAAGCCGAGGGCCAGCGTGCCCGCCGCTGTGTGTCAGCTACGATGACGACGCGGATGAGC
 TCAGCGCTTTCTGCAGCTCCAGGAACCTCACGCGCCTGCCTGATGGAGTCCCGGGCGGCACCCAAGCCCT
 GTGGCTGGACGGCAACAACCTCTCGTCCGTCCTCCCGGCAGCCTTCCAGAACCTCTCCAGCCTGGGCTT
 CTAACCTGCAGGGCGCCAGCTGGGCAGCCTGGAGCCACAGGCGCTGCTGGGCCTAGAGAACCTGTGCC
 ACCTGCACCTGGAGCGGAACAGCTGCGCAGCCTGGCACTCGGCACGTTTGCACACACGCCCGCGCTGGC
 CTGCTCGGCCTCAGCAACAACCGTCTGAGCAGGCTGGAGGACGGGCTCTTCGAGGGCCTCGGCAGCCTC
 TGGGACCTCAACCTCGGCTGGAATAGCCTGGCGGTGCTCCCCGATGCGGGGTTCCGCGGCCTGGGCAGCC
 TGCAGGAGCTGGTGTGGCGGGCAACAGGCTGGCCTACCTGCAGCCCGCGCTCTTCAGCGGCCTGGCCGA
 GCTCCGGGAGCTGGACCTGAGCAGGAACGCGCTGCGGGCCATCAAGGCAACGTGTTTCGTGCAGCTGCC
 CGGCTCCAGAACTCTACCTGGACCGCAACCTCATCGCTGCCGTGGCCCCGGGCGCTTCTGGGCTGA
 AGGCGCTGCGATGGTGGACCTGTCCACAACCGCGTGGCTGGCTCCTGGAGGACAGTTCCCCGGTCT
 GCTGGGCTGCGTGTGCTGCGGCTGTCCACAACGCCATCGCCAGCCTGCGGGCCCGCACCTTCAAGGAC
 CTGCACCTCCTGGAGGAGCTGCAGCTGGCCACAACCGCATCCGGCAGCTGGCTGAGCGCAGCTTTGAGG
 GCCTGGGGCAGCTTGAGGTGCTCAGCTAGACCACAACAGCTCCAGGAGGTCAAGGCGGGCGCTTCTCT
 CGGCCTACCAACGTGGCGGTATGAACCTCTCTGGGAAGTGTCTCCGGAACCTTCCGGAGCAGGTGTTT
 CGGGCCTGGGCAAGCTGCACAGCCTGCACCTGGAGGGCAGCTGCCTGGGACGCATCCGCCCGCACACCT
 TCACCGGCCTCTCGGGGCTCCGCCGACTTTCTCAAGGACAACGGCCTCGTGGGATTGAGGAGCAGAG
 CCTGTGGGGGCTGGCGGAGCTGCTGGAGCTCGACCTGACCTCCAACCAGCTCACGCACCTGCCCCACCGC
 CTCTCCAGGGCCTGGGCAAGCTGGAGTACCTGCTGCTCTCCCGCAACCGCCTGGCAGAGCTGCCGGCGG
 ACGCCCTGGGCCCCCTGCAGCGGCCTTCTGGCTGGACGTCTCGACAACCGCCTGGAGGATTGCCCAA
 CAGCCTTTGGCACCCTGGGGCGGCTGCGCTACCTCAGCCTCAGGAACAACCTACTGCGGACCTTCACG
 CCGCAGCCCCGGCCTGGAGCGCTGTGGCTGGAGGGTAACCCCTGGGACTGTGGCTGCCCTCTCAAGG
 CGCTGCGGGACTTCGCCCTGCAGAACCCAGTGTGTGCCCGCTTCGTCCAGGCCATCTGTGAGGGGGA
 CGATTGCCAGCCGCCCGGTACACCTACAACAACATCACCTGTGCCAGCCGCCCGAGGTCGTGGGGCTC
 GACCTGCGGGACCTCAGCGAGGCCCACTTTGCTCCCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228516 representing NM_001146006
Red=Cloning site Green=Tags(s)


```
MALRKAGDLEPQFTPERFRFCWYQAHSGRALLGPPPQASPPAGGLALALLLSWVALGPRSLEGADPGT
PGEAEGPACPAACVCSYDDDADEL SVFCSRNL TRLPDGVPGGTQALWLDGNLSSVPPAAFQNLSSLGF
LNLQGGQLGSLEPQALLGLENLCHLHLERNQLRSLALGTF AHTPALASLGLSNNRLSRLEDGLFEGLGSL
WDLNLGWNSLAVLPDAAFRGLGSLRELVLAGNRLAYLQPALF SGLAELRELDL SRNALRAIKANVFVQLP
RLQKLYLDRNLI AAVAPGAF LGLKALRWLDL SHNRVAGLLEDTFPGLLGLRVLRLSHNAIASLRPRTFKD
LHFLEELQLGHNRIQLAERSFEGLGQLEVL TLDHNQLQEVKAGAF LGLTNVAVMNL SGNCLRNLPEQVF
RGLGKLSLHLEGSCLGRIRPHTFTGLSGLRRLFLKDNGLVGIEEQSLWGLAELELDLTSNQLTHLPHR
LFQGLGKLEYLLL SRNRLAELPADALGPLQRAF WLDVSHNRLEALPNSLLAPLGRRLRYLSLRNNSLRFTT
PQPPGLERLWLEGNPWDCGCP LKALRDFAL QNPSAVPRFVQA ICEGDDCQPPAYTYNNITCASPEVVGL
DLRDLSEAHFAPC
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8071_c04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001146006

ORF Size: 1929 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146006.2](#)

RefSeq ORF: 1932 bp

Locus ID: 3483

UniProt ID: [P35858](#)


Cytogenetics: 16p13.3

Protein Families: Druggable Genome, Secreted Protein

MW: 70.24 kDa

Gene Summary: The protein encoded by this gene is a serum protein that binds insulin-like growth factors, increasing their half-life and their vascular localization. Production of the encoded protein, which contains twenty leucine-rich repeats, is stimulated by growth hormone. Defects in this gene are a cause of acid-labile subunit deficiency, which manifests itself in a delayed and slow puberty. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RC228516