

## Product datasheet for **RC228522L3V**

### CCDC120 (NM\_001163321) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC120 (NM_001163321) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC120
Synonyms:	JM11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163321
ORF Size:	1983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228522).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163321.1</a> , <a href="#">NP_001156793.1</a>
RefSeq ORF:	2091 bp
Locus ID:	90060
UniProt ID:	<a href="#">Q96HB5</a>
Cytogenetics:	Xp11.23
MW:	70.7 kDa
Gene Summary:	This gene encodes a protein that contains a coiled-coil domain. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. [provided by RefSeq, Nov 2012]


[View online »](#)