

Product datasheet for **RC228524**

ZNF44 (NM_001164276) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF44 (NM_001164276) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF44
Synonyms:	GIOT-2; KOX7; ZNF; ZNF55; ZNF58; ZNF504
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC228524 representing NM_001164276
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCCTGTGTTATGGAACCTTTCTGGGGCTACCCTAAGATGCTGGAAGCTGCCAATCTCATGGAGGGCC
 TAGTGGATATCGGCCCTTGGGTCACTCTCCAGAGGACAGCCTGAGGTGTTGGAGTGGGCCTCCCAAA
 GGATCAGGACTCAGTGGCCTTTGAGGATGTGGCTGTGAACCTCACCCATGAGGAGTGGGCTTTGCTGGGT
 CCATCACAGAAGAATCTCTACAGAGATGTGATGCGAGAAACCATAGGAACCTGAACTGTATAGGAATGA
 AATGGGAAAACGAGAACATTGATGATCAGCACCAAAATCTCAGGAGAAATCCAAGGTGTGATGTGGTAGA
 GAGATTTGGTAAAAGTAAAGATGGTAGTCAGTGTGGAGAAACCTTAAGCCAGATTCGAAATAGTATTGTA
 AACAGAACACTCCCGCCAGAGTAGATGCATGTGGAAGCAGTGTGAATGGAGAAGTCATAATGGGTCATT
 CATCCCTGAATTGCTACATCAGAGTTGATACTGGACACAAACACCGGGAGTGTGATGAATATGCAGAGAA
 GTCATATACACATAAGCAGTGTGGGAAAGGCTTAAGTTATCGCCACTCCTTTCAAACATGTGAAAGGCCT
 CACACTGGAAAGAAACCTATGATTGTAAGGAATGTGGAAAAACCTTCAGTTCTCCTGGAAACCTTCGAA
 GACATATGGTAGTAAAAGGTGGAGATGGACCTTATAAATGTGAATTGTGTGGGAAAGCCTTTTTTGGCC
 CAGTTTATTACGTATGCATGAAAGAACTCACACTGGAGAGAAACCATATGAATGTAAGCAGTGTTCTAA
 GCCTTCCTCTGTTACAGTTCCTATCTAAGACATGAAAAATACACACTGGGGAGAAACCGTATGAATGTA
 AGCAGTGTTCTAAAGCCTTCCTGATTACAGTTCATATCTAAGACATGAAAGAACTCACACTGGAGAGAA
 ACCCTACAAATGTAACAATGTGGGAAAGCCTTCAGTGTTCGGTTCCTTCGAGTACATGAAAGAATT
 CACACTGGAGAGAAACCTATACATGTAACAGTGTGGGAAAGCGTTTGTGATCTTGGAGCTTTCAAA
 GACACATGATAATGCACAGTGGAGATGGACCTCATAAATGTAAGATATGTGGGAAAGCCTTTGATTTTCC
 TGGTTCAGCACGAATTCATGAAGGAACTCACACTCTAGAGAAACCTATGAATGTAAGCAATGTGGGAAA
 TTGTTATCTCATCGCTCAAGCTTTCGAAGACACATGATGGCACACACTGGAGATGGCCCTCATAAATGCA
 CAGTATGTGGGAAAGCCTTTGATTCTCCTAGTGTATTTCAAAGACATGAAAGGACTCACACTGGAGAGAA
 ACCCTATGAATGCAAGCAATGTGGGAAAGCCTTCCTACTTCCAGTTCCTTCGAAACATGAAACAACA
 CACACTGGAGAGCAACCTATAAATGTAATGTGGGAAAGCTTTAGTGATTTATTTTCTTTCAAAGTC
 ATGAAACAACACACAGTGAAGAGGAGCCTTATGAATGTAAGGAGTGTGGGAAAGCATTTAGTTCTTTTAA
 ATACTTTTGTGCGCATGAAAGGACTCACAGTGAAGAAAAATCTTATGAGTGTCAAATTTGTGGCAAAGCC
 TTCAGTCGTTTCAGTTACTTAAAACTCATGAAAGGACTCACACGCAGAGAAGCCATATGAATGTAAGC
 AATGCAGGAAAGCATTCTTTGGCCCTCTTCTTCTAAGACATGAAAGGACTCACACTGGAGAAAGACC
 CTATGAATGTAACACTGTGGTAAAGCCTTCAGTCGTTCCAGTTTCTGTGAGAACATGAAAGAACTCAC
 ACTGGAGAGAAGCCCTATGAATGTAAGGAATGTGGGAAAGCCTTCAGTTCTCTCAGTTCTTTAATAGAC
 ATAAAAGGACACACTGGAAGGATATTCTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC228524 representing NM_001164276
 Red=Cloning site Green=Tags(s)


MALCYGTFWGYPKMLEAANLMEGLVDIGPWVTLPRGQPEVLEWGLPKDQDSVAFEDVAVNFTHEEWALLG
 PSQKNLYRDMRETIRNLNLCIGMKWENQNIDDQHQLRRNPRCDVVERFGKSKDGSQCGETLSQIRNSIV
 NKNTPARVDACGSSVNGEVIMGHSSLNCYIRVDTGHHKRECHYEAKSYTHKQCGKGLSYRHSFQTCERP
 HTGKKPYDCKEKGKTFSSPGNLRHMVVKGGDPYKCELCGKAFFWPSLLRMHERHTHTGEKPYECKQCSK
 AFPVYSSYL RHEKIHTGEKPYECKQCSKAFFDYSSYL RHERHTHTGEKPYCKQCGKAFSVSGSLRVHERI
 HTGEKPYTCKQCGKAFCHLGSFQRHMIMHSGDGPBKCKICGKGDFPGSARIHEGHTHTLEKPYECKQCGK
 LLSHRSSFRHMAHTGDGPHKCTVCGKAFDSPSVFQRHERHTHTGEKPYECKQCGKAFRTSSSLRKHETT
 HTGEQPYKCKCGKAFSDLFSFQSHETTHSEEPYECKEKGKAFSSFKYFCRHERTHSEKSYECQICGKA
 FSRFSYLKTHERTHTAEKPYECKQCRKAFFWPSFLLRHERHTHTGERPYECKHCGKAFSRSSFCREHERTH
 TGEKPYECKEKGKAFSSSFNRHKRTHWKDIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8038_f06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001164276

ORF Size: 1989 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001164276.2
RefSeq ORF:	1992 bp
Locus ID:	51710
UniProt ID:	P15621
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
MW:	76.3 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC228524