

## Product datasheet for **RC228541**

### DISC1 (NM\_001164541) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DISC1 (NM_001164541) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DISC1
Synonyms:	C1orf136; SCZD9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RC228541 representing NM\_001164541  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCCAGGCGGGGTCTCAGGGCGCCACGCCCGCCGGCGGGCGGCGTGTAGCCACCGCGCAGGCA  
GCCGGGATTGCTTACCACCTGCAGCGTGCCTTCGGAGGCGCGGCTGGCACGGAGCCGGCTACATGAG  
AAGCTCGACAGGGCCTGGGATCGGGTCTTTCCCCAGCAGTGGGCACACTGTTCCGGTCCAGGAGGG  
GTGTCTGGCGAGGAGTCCACCCTCGGAGTCCAGGGCCAGACAGTGTGGCCTTGACTCGAGAGGCCTCT  
TGGTCCGGAGCCCTGTTTCCAAGAGTGCAGCAGCCCTACTGTGACCTCTGTGAGAGGAACCTCGGC  
CTTTGGGATTCAGCTCAGAGGTGGCACCAGATTGCCTGACAGGCTTAGCTGGCCGTGTGGCCCTGGGAGT  
GCTGGGTGGCAGCAAGAGTTTGCAGCCATGGATAGTTCTGAGACCCTGGACGCCAGCTGGGAGGCAGCT  
GCAGCGATGGAGCAAGGCGTGTCCGGCCAGCAGGCTCTCTGCCATCAGCAGAGTTGAGTAGCAACAGCTG  
CAGCCCTGGCTGTGGCCCTGAGGTCCCCCAACCCTCCTGGCTCTCACAGTGCCTTTACCTCAAGCTTT  
AGCTTTATTTCGGCTCTCGCTTGGCTCTGCCGGGAACGTGGAGAAGCAGAAGGCTGCCACCATCCAGAG  
AGGCTGAGTCCCATTGCCAGAGCCCCAGGAGATGGGAGCCAAAGCTGCCAGCTTGGACGGGCTCACGA  
GGACCCGCGATGTCTCTCTCGGCCCTTCACTCTTTGGCTACACGGGTCTCTGCAGACTTGGCCAGGCC  
GCAAGGAACAGCTCCAGGCCAGAGCGTGACATGCATTCTTACCAGACATGGACCCTGGCTCCTCCAGTT  
CTCTGGATCCCTCACTGGCTGGCTGTGGTGGTGTGGGAGCAGCGGCTCAGGGGATGCCACTCTTGGGA  
CACCTGTCTCAGGAAATGGGAGCCAGTGTGCGGGACTGCCTGCTGAGAAACCGGAGGCAGATGGAGGTA  
ATATCCTAAGATTAACAACTTACAGAACTTACAGAAAGTGCAGTTGAGAATGATGATTATGATAAAGCTG  
AGACGTTACAACAAAGATTAGAAGACTTGAACAAGAGAAAAATCAGCCTGCACCTTCAACTTCTTCAAG  
GCAGCCAGCTCTTAGCAGTTTCTGGGTACCTGGCAGCACAAAGTCCAGGCTGCCTTGCGCCGTGGGGCC  
ACTCAGCAGGCCAGCGGAGATGACACCCACACCCCACTGAGAATGGAGCCGAGGCTGTTGGAACCCACTG  
CTCAGGACAGCTTGCACGTGTCCATCACGAGACGAGACTGGCTTCTTACGGAAAAGCAGCAGCTACAGAA  
AGAAATCGAAGCTCTCAAGCAAGGATGTTTGTGCTGGAAGCCAAAGATCAACAGCTGAGAAGGGAAATA  
GAGGAGCAAGAGCAGCAACTCCAGTGGCAGGGCTGCGACCTGACCCCACTGGTGGCCAGCTGTCCCTGG  
GTCAGCTGCAGGAGGTCAGCAAGGCTTGCAGGACACCCTGGCCTCAGCCGGTCAAGTCCCTTCCATGC  
AGAGCCACCGAAACCATAAGGAGCCTCCAGGAAAGAATAAAATCCCTCAACTTGTCACTTAAAGAAATC  
ACTACTAAGGTGTGTATGAGTGAGAAATCTGCAGCACCTGAGGAAGAAAGTTAACGATATTGAAACCC  
AACTACCAGCCTTGCTTGAAGCCAAAATGCATGCCATATCAGGAAACCATTTCTGGACGGCTAAAGACCT  
CACCGAGGAGATTAGATCATTAAACATCAGAGAGAGAAGGGCTGGAGGGACTCCTCAGCAAGCTGTTGGTG  
TTGAGTTCCAGGAATGTCAAAAAGCTGGGAAGTGTAAAGAAGATTACAACAGACTGAGAAGAGAAGTGG  
AGCACCAGGAGACTGCCTATGCTGCAAGTGTCCACTGCTTGGGAAAGTGTGGGAAGCTGACTTGGAAAGCT  
TGTCGATTGCTTATCCAGAGCCTACAGCTCCAGGAAGCCAGGGGAAGCCTGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC228541 representing NM\_001164541  
Red=Cloning site Green=Tags(s)


MPGGGPQGAPAAAGGGGVSHRAGSRDCLPPAACFRRRRLARRPGYMRSSSTGPGIGFLSPAVGTLFRFPGG  
 VSGEESHSESRARQCGLDSRGLLVRSVSKSAAAPT VTSVRGTAHFGIQLRGGTRLPDRLSWPCGPGS  
 AGWQQEF AAMDSSSETLDASWEAACSDGARRVRAAGSLPSAELSSNSCSPGCGPEVPPTPPGSHSAFTSSF  
 SFIRLSLGSAGERGEAEGCPPSREAESHCSQSPQEMGAKAASLDGPHEDPRCLSRPFSLLATRVSADLAQA  
 ARNSSRPERDMHSLPDMDPGSSSSLDPSLAGCGDGGSSGSDAHSWDTLLRKWEPVLRDCLLRNRRQMEV  
 ISLRRLKQLKLEDVENDDYDKAETLQQRLLEDLEQEKISLHFQLPSRQPALSSFLGHLAAQVQAALRRGA  
 TQQASGDDTHTPLRMEPRLLEPTAQDSLHVSITRRDWLLQEKKLQKEIEALQARMFVLEAKDQQLRREI  
 EEQEQLQWQGC DL TPLVGQLSLGQLQEVSKALQDTLASAGQIPFHAEPPE TIRSLQERIKSLNLSLKEI  
 TTKVCMSEKFCSTLRKKVNDIETQLPALLEAKMHAISGNHFWTAKDL TEEIRSLT SEREGLEGLLSKLLV  
 LSSRNVKKLSVKEDYNRLRREVEHQETAYAASVHCLGKCGKLTWKLVDCLSRAYSRRKPGEACL

TRTRPLEQKLISEEDLANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8073\\_d01.zip](https://cdn.origene.com/chromatograms/mk8073_d01.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**




**ACCN:** NM\_001164541

**ORF Size:** 2085 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001164541.2</u>
<b>RefSeq ORF:</b>	2088 bp
<b>Locus ID:</b>	27185
<b>UniProt ID:</b>	<u>Q9NRI5</u>
<b>Cytogenetics:</b>	1q42.2
<b>MW:</b>	75.4 kDa
<b>Gene Summary:</b>	This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

**Product images:**


Circular map for RC228541