

Product datasheet for **RC228571**

PDE4D (NM_001165899) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE4D (NM_001165899) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PDE4D
Synonyms:	ACRDYS2; DPDE3; HSPDE4D; PDE4DN2; PDE43; STRK1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC228571 representing NM_001165899
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAAGAAATACCTGTGATTTGCTTTCTCGGAGCAAAAGTGCCTCTGAGGAAACTACATTCCAGTA
 ATGAAGAGGAAGACCCTTCCGCGGAATGGAACCTATCTTGCCGGAGACTTTCATGTCGCAATATTCA
 GCTTCCCCTCTCGCCTTCAAGACAGTTGGAACAAGCTGACTTGAAAAGTGAATCAGAGAATTTCAACGA
 CCAACCAGCCTCCCCTGAAGATTCTGCCGCTGATTGCTACTTCTGCAGAATCCAGTGGTTTTGATG
 TGGACAATGGCACATCTGCGGGACGGAGTCCCTGGATCCATGACCAGCCAGGATCCGGGCTAATTCT
 CCAAGCAAATTTGTCCACAGTCAACGACGGGAGTCCCTTCTGTATCGATCCGACAGCGATTATGACCTC
 TCTCAAAGTCTATGTCGGAACTCCTCCATTGCCAGTGATATACACGGAGATGACTTGATTGTGACTC
 CATTGCTCAGGCTTGCCAGTCTGCGAACTGTACGAAACAACCTTGTGCATTAATAATTTGCAAGA
 TCGAGCACCTAGCAAAAGATCACCCATGTGCAACCAACCATCCATCAACAAAGCCACCATAACAGAGGAG
 GCCTACCAGAAACTGGCCAGCGAGACCCTGGAGGAGCTGGACTGGTGTCTGGACCAGCTAGAGACCCTAC
 AGACCAGGCACTCCGTCAGTGAGATGGCCTCCAACAAGTTTAAAAGGATGCTTAATCGGGAGCTCACCCA
 TCTCTGAAATGAGTCCGTCTGAAATCAAGTGTGAGAGTTTATATCAAACACATTTAGATAAGCAA
 CATGAAGTGGAAATTCCTTCTCCAACCTCAGAAGGAAAAGGAGAAAAAGAAAGACCAATGTCTCAGATCA
 GTGGAGTCAAGAAATGATGCACAGCTCTAGTCTGACTAATCAAGTATCCCAAGGTTTGGAGTTAAAAC
 TGAACAAGAAGATGTCCTTGCCAAGGAAGTGAAGATGTGAACAAATGGGGTCTTCATGTTTTGAGAATA
 GCAGAGTTGTCTGGTAACCGGCCCTTACTGTTATCATGCACACCATTTTCAGGAACGGGATTTATTA
 AAACATTTAAAATTCAGTAGATACTTTAATTACATATCTTATGACTCTCGAAGACCATTACCATGTCTGA
 TGTGGCCTATCAACAATATCCATGCTGCAGATGTTGTCCAGTCTACTCATGTGCTATTATCTACACCT
 GCTTTGGAGGCTGTGTTACAGATTTGGAGATTCTGCAGCAATTTTGGCAGTGAATACATGATGTAG
 ATCATCTGTTGTGTCATCAATTTCTGATCAATACAACTCTGAACTTGCCTTGATGTACAATGATTC
 CTCAGTCTTAGAGAACCATCATTGGCTGTGGCTTTAAATTGCTTCAGGAAGAAAAGTGTGACATTTTC
 CAGAAATTTGACCAAAAAACAAAGACAATCTTAAAGGAAAATGGTCATTGACATCGTACTTGAACAGATA
 TGTCAAAACACATGAATCTACTGGCTGATTTGAAGACTATGGTTGAAACTAAGAAAGTGAACAGCTCTGG
 AGTTCTTCTTCTGATAATTATCCGATAGGATTCAGGTTCTTCAAGATATGGTGCAGTGTGCAGATCTG
 AGCAACCAACAAGCCTCTCCAGCTGTACCGCCAGTGGACGGACCGGATAATGGAGGATTTCTCCGCC
 AAGGAGACCGAGAGAGGGAACGTGGCATGGAGATAAGCCCATGTGTGACAAGCACAATGCTTCCGTGGA
 AAAATCACAGGTGGGCTTCATAGACTATATTGTTTCATCCCTCTGGGAGACATGGGCAGACCTCGTCCAC
 CCTGACGCCAGGATATTTGGACACTTTGGAGGACAATCGTGAATGGTACCAGAGCACAATCCCTCAGA
 GCCCTCTCTGCACCTGATGACCCAGAGGAGGGCCGGCAGGGTCAAACCTGAGAAATTCAGTTTGAAC
 AACTTTAGAGGAAGATGGTGTGAGTCAACACGAAAAGGACAGTGGCAGTCAAGTGGAAAGAAGACTAGC
 TGCAGTACTCCAAGACTCTTTGTACTCAAGACTCAGAGTCTACTGAAATTTCCCTTGTGTAACAGGTTG
 AAGAGGAGGCAGTAGGGGAAGAAGAGGAAAAGCCAGCCTGAAGCCTGTGTATAGATGATCGTTCTCTGTA
 CACG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001165899.2
RefSeq ORF:	2247 bp
Locus ID:	5144
UniProt ID:	Q08499
Cytogenetics:	5q11.2-q12.1
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism
MW:	84.5 kDa
Gene Summary:	This gene encodes one of four mammalian counterparts to the fruit fly 'dunce' gene. The encoded protein has 3',5'-cyclic-AMP phosphodiesterase activity and degrades cAMP, which acts as a signal transduction molecule in multiple cell types. This gene uses different promoters to generate multiple alternatively spliced transcript variants that encode functional proteins.[provided by RefSeq, Sep 2009]

Product images:



Circular map for RC228571