

## Product datasheet for **RC228576L4V**

### APBB2 (NM\_001166050) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	APBB2 (NM_001166050) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APBB2
Synonyms:	FE65L; FE65L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166050
ORF Size:	2274 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC228576).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001166050.1</a>
RefSeq ORF:	2277 bp
Locus ID:	323
UniProt ID:	<a href="#">Q92870</a>
Cytogenetics:	4p14-p13
Protein Families:	Transcription Factors
MW:	83.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene interacts with the cytoplasmic domains of amyloid beta (A4) precursor protein and amyloid beta (A4) precursor-like protein 2. This protein contains two phosphotyrosine binding (PTB) domains, which are thought to function in signal transduction. Polymorphisms in this gene have been associated with Alzheimer's disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]