

Product datasheet for **RC228598**

CTNNA2 (NM_001164883) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTNNA2 (NM_001164883) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CTNNA2
Synonyms:	CAP-R; CAPR; CDCBM9; CT114; CTNR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC228598 representing NM_001164883
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACTTCGGCAACTTCACCTATCATTCTGAAATGGGACCCAAAAGTTTGGAAATCCGGACGCTAACAG
 TGGAAAGGCTGTTGGAGCCACTTGTACACAGGTGACTACACTTGTCAACACAAGCAACAAAGGCCCATC
 TGGTAAAAAGAAAGGGAGGTCAAAGAAAGCCCATGTACTAGCTGCCTCTGTAGAGCAAGCCACTCAGAAT
 TTCCTGAAAAAGGGTGAACAGATCGCTAAGGAGAGTCAAGATCTCAAAGAAGAGTTGGTGCTGTGTAG
 AGGATGTGCGCAAACAAGGTGAGACGATGCGGATGCGCTCCTCGAGTTTGCAGATGACCCTTGCTCGTC
 GGTAAAGCGCGGCACCATGGTACGGGCGCAAGGGCTTTGCTCTCCGCGGTGACACGCTTACTCATCCTG
 GCGGACATGGCAGATGTCATGAGACTTTTATCCCATCTGAAAATTGTGGAAGAGGCCCTGGAAGCTGTCA
 AAAATGCTACAAATGAGCAAGACCTTGCAAACCGTTTTAAAGAGTTTGGGAAAGAGATGGTAAACTTAA
 CTATGTAGCAGCAAGAAGACAACAGGAGCTGAAGGATCCTCACTGTCGGGATGAGATGGCAGCCGCCGA
 GGGGCTCTGAAGAAGAATGCCACAATGCTGTACACGGCCTCTCAAGCATTCTCCGCCACCCAGATGTCCG
 CCGCTACGAGAGCCAACCGAGATTATGTGTTCAAACAAGTCCAGGAGGCCATCGCCGGCATCTCCAATGC
 TGCTCAAGCTACCTCGCCACTGACGAAGCCAAGGGCCACACGGGCATCGGCGAGCTGGCTGCGGCTCTT
 AATGAGTTTGACAATAAGATTATCCTGGACCCCATGACGTTTACGCGAGGCCAGGTTCCGGCCGCTCCCTGG
 AGGAGAGGCTGGAGAGCATCATCAGCGGCGCAGCGCTGATGGCCGACTCCTCCTGCACGCGAGACGACCG
 GCGCGAGAGGATCGTGGCGGAGTGAACGCCGTGCGGCAGGCGCTCCAGGACCTGCTCAGCGAGTACATG
 AATAACTGGAAGGAAAGAAAAAGGAGATCCTCTCAACATTGCGATTGATAAGATGACTAAGAAAAACA
 GAGATCTAAGGAGACAGCTTCGAAAGCAGTGATGGATCACATATCTGACTCTTTCTGGAACCAATGTT
 TCCTTTGCTAGTTCTCATTGAGGCTGCAAAGAGCGGAAATGAAAAGGAAGTGAAGAATATGCCCAAGTT
 TTCGTGAGCATGCCAACAACTGGTAGAGTTGCCAATTTGGCCTGTTCCATCTCCAACAATGAAGAAG
 GGGTGAATTAGTTCCGATGGCAGCCACCAGATTGACAGCCTGTGTCCCAGGTCATCAATGCCGCTCT
 GACACTGGCTGCCCGCCACAGAGCAAAGTTGCTCAGGATAACATGGACGCTTTCAAAGACCAGTGGGAG
 AAGCAGGTCAGTGTGACAGAGGCCGTGGATGACATCACCTCAGTGGATGACTTCTCTGTCTCAG
 AAAATCACATCTGGAGGATGTGAACAAGTGTGTGATAGCCCTCAAGAGGGCGATGTGGACACTCTGGA
 CCGGACTGCAGGGGCCATCAGGGGCCGGCAGCTCGAGTCATACATCATCAATGCTGAGATGGAGAAC
 TATGAAGCTGGGTTTATACTGAGAAGGTGTTGGAAGCTACAAAATTGCTTTCTGAAACAGTGATGCCAC
 GCTTCGCTGAACAAGTAGAGTTGCCATTGAAGCCCTGAGTGCCAACGTTCTCAACCGTTTGAGGAGAA
 TGAGTTCATCGATGCCTCTCGCCTGGTGTATGATGGCGTTCGGGACATCAGAAAGGCTGTGCTGATGATC
 AGGACCCAGAAAGAACTAGAGGATGATTCTGACTTTGAGCAGGAAGATTATGATGTGCGTAGCAGGACAA
 GTGTTCCAGACTGAGGATGACCAGCTCATTGAGGGCAGAGCGCACGGGCCATCATGGCGCAACTACCGCA
 GGAGGAGAAGGCAAAAATAGCTGAGCAGGTGGAGATATCCATCAAGAGAAAAGCAAGCTGGATGCAGAA
 GTGGCCAAATGGGACGACAGCGGCAATGATATCATTGACTGGCCAAGCAGATGTGTATGATCATGATGG
 AAATGACAGACTTCACAAGAGGCAAAGGCCATTGAAAAATACATCTGATGTCTAATGCTGCCAAGAA
 AATTGCCGAAGCAGGTTCTCGAATGGACAAATTAGCTCGTGTGGCTGATCAGCTGGACAGTGCCACA
 TCGCTTATCCAGGCAGCTAAAAACCTGATGAATGCTGTTGTCTCACGGTGAAGCATCCTATGTGGCCT
 CAACCAAATACGAAAGTCTATGGGACAGCAGCTGTCAACTCACCTGTTGTGCTTGGAAAGATGAAGGC
 TCCAGAGAAGAAGCCCTTGTGAAGAGAGAAAAGCCTGAAGAATTCCAGACACGAGTTCGACGAGGTTCT
 CAGAAGAAACACATTTCCGCTGTACAGGCTTAAAGTGAATTCAAAGCAATGGATTCTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

ACCN:	NM_001164883
ORF Size:	2580 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001164883.1 , NP_001158355.1
RefSeq ORF:	2583 bp
Locus ID:	1496
UniProt ID:	P26232
Cytogenetics:	2p12
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Endometrial cancer, Leukocyte transendothelial migration, Pathways in cancer, Tight junction
MW:	95.4 kDa
Gene Summary:	May function as a linker between cadherin adhesion receptors and the cytoskeleton to regulate cell-cell adhesion and differentiation in the nervous system (By similarity). Required for proper regulation of cortical neuronal migration and neurite growth (PubMed:30013181). It acts as negative regulator of Arp2/3 complex activity and Arp2/3-mediated actin polymerization (PubMed:30013181). It thereby suppresses excessive actin branching which would impair neurite growth and stability (PubMed:30013181). Regulates morphological plasticity of synapses and cerebellar and hippocampal lamination during development. Functions in the control of startle modulation (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC228598